

Product datasheet for **TP311896L**

CCDC43 (NM_144609) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 43 (CCDC43), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211896 representing NM_144609 Red =Cloning site Green =Tags(s)

MAAPSEVAAIAPGEGDGGGGFGSWLDGRLEALGVDRVYGAYILGILQEEEEEEKLDALQGILSAFLEE
DSLLNICKEIVERWSETQNVVTKVKKKEDEVQAIATLIEKQAQIVVKPRMVSEEEKQRKAALLAQYADVTD
EEDEADEKDDSGATTMNIGSDKLLFRNTNVEDVLNARKLERDSLDESQRKKEQDKLQRRERDKLAKQERK
EKEKKRTQRGERKR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	25.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_653210</u>
Locus ID:	124808



[View online »](#)

UniProt ID: [Q96MW1](#)
RefSeq Size: 2082
Cytogenetics: 17q21.31
RefSeq ORF: 672

Product images:



Coomassie blue staining of purified CCDC43 protein (Cat# [TP311896]). The protein was produced from HEK293T cells transfected with CCDC43 cDNA clone (Cat# [RC211896]) using MegaTran 2.0 (Cat# [TT210002]).