

## Product datasheet for **TP311900L**

### INPP4A (NM\_004027) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins  
**Description:** Recombinant protein of human inositol polyphosphate-4-phosphatase, type I, 107kDa (INPP4A), transcript variant a, 1 mg  
**Species:** Human  
**Expression Host:** HEK293T  
**Expression cDNA Clone or AA Sequence:** >RC211900 representing NM\_004027  
**Red**=Cloning site **Green**=Tags(s)

MTAREHSPRHGARARAMQRASTIDVAADMLGLSLAGNIQDPDEPILEFSLACSELHTPSLDRKPNSFVAV  
SVTTPPQAFWTKHAQTEIIEGTNNPIFLSSIAFFQDSLINQMTQVKLSVYDVKDRSQGTMVLLGSGTFIV  
KDLLQDRHRLHLTLRSAESDRVGNITVIGWQMEEKSDQRPPVTRSVDTVNGRMVLPVDESALTEALGIRS  
KYASLRKDTLLKSVFGAICRMYRFPTTDGNHLRILEQMAESVLSLHVPRQFVKLLLEEDAARVCELEEL  
GELSPCWESLRRQIVTYQTIILTYQENLTDLHQYRGPFSFKASSLKADKKLEFVPTNLHIQRMVQDDGG  
SDQNYDIVTIGAPAAHCQGFKSGGLRKKLHKFEETKKHFEECCTSSGCQSIYIPQDVRAKEIIAQINT  
LKTQVSYAERLSRAAKDRSATGLERTLAILADKTRQLVTVCDCKLLANSIHGLNAARPDYIASKASPTS  
TEEEQVMLRNDQDTLMARWTGRNSRSSLQVDWHEEEWEKVWLVNVDKSLECIQRVDKLLQKERLHGECE  
DVFP CAGSCTSKKDCSPPEESSPGEWSEALYPLLTLDVCVAMMSDKAKKAMVFLMQDSAPTIAATYLS  
LQYRRDVVFCQTLTALICGFIIKLRNCLHDDGFLRQLYTIGLLAQFESLLSTYGEELAMLEDMSLGIMDL  
RNVTFKVTQATSSASADMLPVITGNRDGFNVRVPLPGPLFDALPREIQSGMLLRVQPVLFVNGINEQRTL  
AERFGDTSLQEVINVESLVRNLSYFEQFKEVLPEDCLPRSRSTCLPELLRFLGQNVHARKNKNVDILWQ  
AAEICRRLNGVRFTSCKSAKDRTAMSVTLEQCLILQHEHGMAPQVFTQALECMRSEGCRRRENTMKNVGSR  
KYAFNSLQLKAFPKHYRPPPEGTYGKVET

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


**Tag:** C-Myc/DDK  
**Predicted MW:** 105.5 kDa  
**Concentration:** >0.05 µg/µL as determined by microplate BCA method  
**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining  
**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_004018</a>
<b>Locus ID:</b>	3631
<b>UniProt ID:</b>	<a href="#">Q96PE3</a>
<b>RefSeq Size:</b>	3217
<b>Cytogenetics:</b>	2q11.2
<b>RefSeq ORF:</b>	2814
<b>Synonyms:</b>	INPP4; TVAS1
<b>Summary:</b>	This gene encodes an Mg <sup>++</sup> independent enzyme that hydrolyzes the 4-position phosphate from the inositol ring of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate, and inositol 3,4-bisphosphate. Multiple transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Aug 2008]
<b>Protein Families:</b>	Transmembrane
<b>Protein Pathways:</b>	Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system

### Product images:



Coomassie blue staining of purified INPP4A protein (Cat# [TP311900]). The protein was produced from HEK293T cells transfected with INPP4A cDNA clone (Cat# [RC211900]) using MegaTran 2.0 (Cat# [TT210002]).