

Product datasheet for TP312051

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PAX2 (NM_003987) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human paired box 2 (PAX2), transcript variant a, 20 μg

Species: Human Expression Host: HEK293T

Expression cDNA >RC212051 representing NM_003987 Clone or AA Sequence: Red=Cloning site Green=Tags(s)

MDMHCKADPFSAMHPGHGGVNQLGGVFVNGRPLPDVVRQRIVELAHQGVRPCDISRQLRVSHGCVSKILG RYYETGSIKPGVIGGSKPKVATPKVVDKIAEYKRQNPTMFAWEIRDRLLAEGICDNDTVPSVSSINRIIR TKVQQPFHPTPDGAGTGVTAPGHTIVPSTASPPVSSASNDPVGSYSINGILGIPRSNGEKRKRDEVEVYT DPAHIRGGGGLHLVWTLRDVSEGSVPNGDSQSGVDSLRKHLRADTFTQQQLEALDRVFERPSYPDVFQAS EHIKSEQGNEYSLPALTPGLDEVKSSLSASTNPELGSNVSGTQTYPVVTGRDMASTTLPGYPPHVPPTGQ GSYPTSTLAGMVPGSEFSGNPYSHPQYTAYNEAWRFSNPALLSSPYYYSAAPRSAPAAAAAAYDRH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 44.5 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 003978

Locus ID: 5076





UniProt ID: Q02962

RefSeq Size: 4276

Cytogenetics: 10q24.31 RefSeq ORF: 1248

Synonyms: FSGS7; PAPRS


Summary: PAX2 encodes paired box gene 2, one of many human homologues of the Drosophila

melanogaster gene prd. The central feature of this transcription factor gene family is the conserved DNA-binding paired box domain. PAX2 is believed to be a target of transcriptional supression by the tumor suppressor gene WT1. Mutations within PAX2 have been shown to result in optic nerve colobomas and renal hypoplasia. Alternative splicing of this gene results in

multiple transcript variants. [provided by RefSeq, Dec 2014]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified PAX2 protein (Cat# TP312051). The protein was produced from HEK293T cells transfected with PAX2 cDNA clone (Cat# [RC212051]) using MegaTran 2.0 (Cat# [TT210002]).