

## Product datasheet for **TP312051M**

### **PAX2 (NM\_003987) Human Recombinant Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human paired box 2 (PAX2), transcript variant a, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC212051 representing NM_003987
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MDMHCKADPFPSAMHPGHGGVNLGGVFNVRPLPDVWRQRIVELAHQGVPCDISRQLRVSHGCVSKILG  
RYYETGSIKPGVIGGSKPKVATPKVVDKIAEYKRQNPTMFAWEIRDLLAEGICDNDTVPSVSSINRIIR  
TKVQQPFHPTPDGAGTGVTAPGHTIVPSTASPPVSSASNDPVGSYSINGILGIPRSNGEKRRKRDVEVYT  
DPAHIRGGGLHLVWTLRDVSEGSVPNGDSQSGVDSLRLKHLRADTFTQQLEALDRVFERPSYDPVFAQS  
EHIKSEQGNEYSLPALTPGLDEVKSSLSASTNPELGSNVSGTQTPVVTGRDMASSTLPGYPPHPVPTGQ  
GSYPTSTLAGMVPVGFSEFSGNPNYSHPQYTAYNEAWRFSNPALLSSPYYSAPRSAPAAAAAYDRH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	44.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_003978</a>
Locus ID:	5076



[View online »](#)

UniProt ID: [Q02962](#)

RefSeq Size: 4276

Cytogenetics: 10q24.31


RefSeq ORF: 1248

Synonyms: FSGS7; PAPRS

**Summary:** PAX2 encodes paired box gene 2, one of many human homologues of the *Drosophila melanogaster* gene *prd*. The central feature of this transcription factor gene family is the conserved DNA-binding paired box domain. PAX2 is believed to be a target of transcriptional suppression by the tumor suppressor gene WT1. Mutations within PAX2 have been shown to result in optic nerve colobomas and renal hypoplasia. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2014]

**Protein Families:** Druggable Genome

### Product images:



Coomassie blue staining of purified PAX2 protein (Cat# [TP312051]). The protein was produced from HEK293T cells transfected with PAX2 cDNA clone (Cat# [RC212051]) using MegaTran 2.0 (Cat# [TT210002]).