



UniProt ID: [Q9NY12](#)

RefSeq Size: 1021

Cytogenetics: 4q25


RefSeq ORF: 651

Synonyms: NOLA1

**Summary:** This gene is a member of the H/ACA snoRNPs (small nucleolar ribonucleoproteins) gene family. snoRNPs are involved in various aspects of rRNA processing and modification and have been classified into two families: C/D and H/ACA. The H/ACA snoRNPs also include the DKC1, NOLA2 and NOLA3 proteins. These four H/ACA snoRNP proteins localize to the dense fibrillar components of nucleoli and to coiled (Cajal) bodies in the nucleus. Both 18S rRNA production and rRNA pseudouridylation are impaired if any one of the four proteins is depleted. These four H/ACA snoRNP proteins are also components of the telomerase complex. The encoded protein of this gene contains two glycine- and arginine-rich domains and is related to *Saccharomyces cerevisiae* Gar1p. Two splice variants have been found for this gene. [provided by RefSeq, Jul 2008]

**Protein Families:** Stem cell - Pluripotency

### Product images:



Coomassie blue staining of purified GAR1 protein (Cat# [TP312074]). The protein was produced from HEK293T cells transfected with GAR1 cDNA clone (Cat# [RC212074]) using MegaTran 2.0 (Cat# [TT210002]).