

Product datasheet for **TP312119L**

CDCP2 (NM_201546) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CUB domain containing protein 2 (CDCP2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212119 representing NM_201546 Red =Cloning site Green =Tags(s)

MLAEWGACLLLAVALLGPGQLQAQAMEGVKCGGVLSAPSGNFSSPNFPRLYPYNTECSWLIVVAEGSSVLL
TFHAFDLEYHDTCSFDFLEIYNGASPDKGNLLGRFCGKVPPPPFTSSWHVMSVIFHSDKHVASHGFSAGY
QKDVCGGVLTGLSGVLTSPPEYNNYPNSMECHWWIRAAGPAHVKLVFVDFQVEGNEECTYDYVAVLGGP
G

PTRGHHYCGSTRPPTLVSLGHELQVVFKSDFNIGGRGFKAYYFSGECQEVYMMAMRGNFSSPQYPSSYPNN
IRCHWTIRLPPGYQVKVFFLDLDLEPNLSLTKTCDFDHLAAFDGAEEAPLLGNWCGHHLPPPVTSSHNQ
LLLLLHTDRSTTRRGFSVAYIGGQLGCGSGSTEGERGEALQPQSLQSPSSIPPVCPAPPMNGLLQLLLHWL
HPCPLSGPLRLDGTAPACFHYCRASFPSF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	48.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: NP_963840
 Locus ID: 200008
 UniProt ID: Q5VXM1
 RefSeq Size: 2723
 Cytogenetics: 1p32.3
 RefSeq ORF: 1347
 Protein Families: Druggable Genome

Product images:


Coomassie blue staining of purified CDCP2 protein (Cat# [TP312119]). The protein was produced from HEK293T cells transfected with CDCP2 cDNA clone (Cat# [RC212119]) using MegaTran 2.0 (Cat# [TT210002]).