

## Product datasheet for TP312143

### Estrogen Related Receptor gamma (ESRRG) (NM\_206594) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human estrogen-related receptor gamma (ESRRG), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212143 representing NM_206594 Red=Cloning site Green=Tags(s)

MSNKDRHIDSSCSSFIKTEPSSPASLTDSVNHHSPPGGSSDASGSYSSTMNGHQNGLDSPPLYPSAPILGG  
SGPVRKLYDDCSSTIVEDPQTKCEYMLNSMPKRLCLVCGDIASGYHYGVASCEACKAFFKRTIQGNIEYS  
CPATNECEITKRRRKSCQACRFMKCLKVGMLEKEGVRLDRVRGGRQKYKRRIDAENSPYLNPLVQPAKKP  
YNKIVSHLLVAEPEKIYAMPDPTVPDSDIKALTTLCDLADRELVIIGWAKHIPGFSTLSLADQMSLLQS  
AWMEILILGWYRSLSFEDLVYADDYIMDEDQSKLAGLLDLNNAILQLVKKYKSMKLEKEEFVTLKAIA  
LANSDSMHIEDVEAVQKLQDVLHEALQDYEAGQHMEDPRRAGKMLMLTPLLRLQTSTKAVQHFYNIKLEGK  
VPMHKLFLMLEAKV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	48.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	WB positive control (PMID: <a href="#">25496115</a> ) WB positive control (PMID: <a href="#">27363015</a> )
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.




[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_996317</a>
<b>Locus ID:</b>	2104
<b>UniProt ID:</b>	<a href="#">P62508</a> , <a href="#">F1D8R6</a>
<b>RefSeq Size:</b>	5369
<b>Cytogenetics:</b>	1q41
<b>RefSeq ORF:</b>	1305
<b>Synonyms:</b>	ERR-gamma; ERR3; ERRg; ERRgamma; NR3B3
<b>Summary:</b>	<p>This gene encodes a member of the estrogen receptor-related receptor (ESRR) family, which belongs to the nuclear hormone receptor superfamily. All members of the ESRR family share an almost identical DNA binding domain, which is composed of two C4-type zinc finger motifs. The ESRR members are orphan nuclear receptors; they bind to the estrogen response element and steroidogenic factor 1 response element, and activate genes controlled by both response elements in the absence of any ligands. The ESRR family is closely related to the estrogen receptor (ER) family. They share target genes, co-regulators and promoters, and by targeting the same set of genes, the ESRRs seem to interfere with the ER-mediated estrogen response in various ways. It has been reported that the family member encoded by this gene functions as a transcriptional activator of DNA cytosine-5-methyltransferases 1 (Dnmt1) expression by direct binding to its response elements in the DNMT1 promoters, modulates cell proliferation and estrogen signaling in breast cancer, and negatively regulates bone morphogenetic protein 2-induced osteoblast differentiation and bone formation. Multiple alternatively spliced transcript variants have been identified, which mainly differ at the 5' end and some of which encode protein isoforms differing in the N-terminal region. [provided by RefSeq, Aug 2011]</p>
<b>Protein Families:</b>	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors

Product images:



Western blot analysis of basal estrogen-related receptor gamma (ERRgamma) expression in different mammary epithelial and breast cancer cell lines. OriGene TP312143 serves as a positive control for ERRgamma. Figure cited from Oncotarget, PMID: 27363015



Coomassie blue staining of purified ESRRG protein (Cat# TP312143). The protein was produced from HEK293T cells transfected with ESRRG cDNA clone (Cat# [RC212143]) using MegaTran 2.0 (Cat# [TT210002]).