

## Product datasheet for **TP312193M**

### PTP4A3 (NM\_032611) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein tyrosine phosphatase type IVA, member 3 (PTP4A3), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212193 representing NM_032611 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MARMNRPAPVEVSYKHMRLITHNPTNATLSTFIEDLKKYGATTVVRVCEVTDKTPLEKDGITVVDWPF DDGAPPPGKVVEDWLSLVKAKFCEAPGSCVAVHCVAGLGRAPVLVALALIESGMKYEDIAIQFIRQKRRGA INSKQLTYLEKYRPKQLRFKDPHTHKTRCCVM  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	19.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_116000</a>
Locus ID:	11156
UniProt ID:	<a href="#">O75365</a>



[View online »](#)

RefSeq Size: 1396

Cytogenetics: 8q24.3


RefSeq ORF: 519

Synonyms: PRL-3; PRL-R; PRL3

**Summary:** This gene encodes a member of the protein-tyrosine phosphatase family. Protein tyrosine phosphatases are cell signaling molecules that play regulatory roles in a variety of cellular processes. Studies of this class of protein tyrosine phosphatase in mice demonstrates that they are prenylated in vivo, suggesting their association with cell plasma membrane. The encoded protein may enhance cell proliferation, and overexpression of this gene has been implicated in tumor metastasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]

**Protein Families:** Druggable Genome, Phosphatase

### Product images:



Coomassie blue staining of purified PTP4A3 protein (Cat# [TP312193]). The protein was produced from HEK293T cells transfected with PTP4A3 cDNA clone (Cat# [RC212193]) using MegaTran 2.0 (Cat# [TT210002]).