

Product datasheet for **TP312200M**

C19orf44 (NM_032207) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 19 open reading frame 44 (C19orf44), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212200 representing NM_032207 Red =Cloning site Green =Tags(s)

MASARKASRPMRDVFGDFSDVSLEDSTMEIRNFQISRNLTKIAPGHSRFLKRNQTLDEKHLLKLENPVL
GSGPRLASCRPPTTASRIRANAALMKLAQLETRIMNRKLQRNLSDTESDMSMTADAGLPKRADRLSGGAL
ELASQNTDKTSQNQARELPVTENNAQNAKVSRLKKKQAPVENISPEAPAGKERTLQTPKQKEPARTFDS
PDSDEEEMKVLLGSLMDSSREKNTNQGFSSANVSEEEERKLFSVPSQLRAFTVPSVELSSAKPSQTSHP
TSLAADRTLHSTRSRADYPQSHVSSDTASHTPSVSITGAFSNSVSLKMGHVKLVSPPGRSEAEVDEPVS
EGADDSLDEFINILSLDGLAPAVSENSDLEQEEESAQRQKTAGKIFRAEASTGQDAPRQAQARSWASQG
KAASAEGDESEVSEHLSASSASAIQDSTSSMQPPSEAPMVNTVSSAYSEDFENSPSLTASEPTAHSKES
LDRTLDALSESSSVKTDLPQTAESRKSGRHVTRVLVKDTAVQTPDPAFTYEWTKVASMAAMGPALGGA
YVDPTPIANHVISADAIEALTAYSPAVLALHDVLKQQLSLTQQFIQASRHLHASLLRSLDADSFHYHTLE
EAKEYIRCHRPAPLTMEDALEEVNKL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	71.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_115583</u>
Locus ID:	84167
UniProt ID:	<u>Q9H6X5</u>
RefSeq Size:	3345
Cytogenetics:	19p13.11
RefSeq ORF:	1971

Product images:


Coomassie blue staining of purified C19orf44 protein (Cat# [TP312200]). The protein was produced from HEK293T cells transfected with C19orf44 cDNA clone (Cat# [RC212200]) using MegaTran 2.0 (Cat# [TT210002]).