

Product datasheet for **TP312208**

CNN2 (NM_201277) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens calponin 2 (CNN2), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC212208 representing NM_201277
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MSSTQFNKGPSYGLSAEVKNRLLSKYDPQKEAELRTWIEGLTGLSIGPDFQKGLKDGTLCTLMNKLQPG
SVPKINRSMQNWHQLENLSNFIKAMVSYGMNPVDLFEANDLFESGNMTQVQVSLLAGKMGTKNCASQS
GMTAYGTRRHLYDPKNHILPPMDHSTISLQMGTKNCASQVGMTAPGTRRHLYDTKLGTDKCDNSSMSLQM
GYTQGANQSGQVFLGRQIYDPKYCPQGTVADGAPSGTGDCPDPGEVPEYPPYYQEEAGY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	29.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_958434
Locus ID:	1265
UniProt ID:	Q99439



[View online »](#)


RefSeq Size: 2361

Cytogenetics: 19p13.3

RefSeq ORF: 810

Summary: The protein encoded by this gene, which can bind actin, calmodulin, troponin C, and tropomyosin, may function in the structural organization of actin filaments. The encoded protein could play a role in smooth muscle contraction and cell adhesion. Several pseudogenes of this gene have been identified, and are present on chromosomes 1, 2, 3, 6, 9, 11, 13, 15, 16, 21 and 22. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jan 2015]

Product images:



Coomassie blue staining of purified CNN2 protein (Cat# TP312208). The protein was produced from HEK293T cells transfected with CNN2 cDNA clone (Cat# [RC212208]) using MegaTran 2.0 (Cat# [TT210002]).