

Product datasheet for TP312231L

SMARCA6 (HELLS) (NM_018063) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human helicase, lymphoid-specific (HELLS), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212231 representing NM_018063 Red=Cloning site Green=Tags(s)

MPAERPAGSGGSEAPAMVEQLDTAVITPAMLEEEEQLEAAGLERERKMLEKARMSWDRESTEIRYRRLQH
LLEKSNIYSKFLLTKMEQQLEEQKKKEKLERKKESLKVKKGKNSIDASEEKPMRKKRGREDESYNISE
VMSKEEILSVAKKNKKENEDENSSSTNLCVEDLQKNKDSNSIIKDRLSETRQNTKFFFDPPVRKCNQGPV
PFQQPKHFTGGVMRWYQVEGMEWLRMLWENGINEILADEMGLGKTVQCIATIALMIQRGVPGPFLVCGPL
STLPNWMAEFKRFTPDIPITMLYHGTQEERQKLVRNIIYKRKGTLQIHPVVITSFEIAMRDRNALQHCYWKY
LIVDEGHRIKNMKRLIRELKRFNADNKLLLTGTPLQNNLSELWSLLNLLPDVDFDLKSFESWFDITSL
SETAEDIIAKEREQNVLHMLHQILTPFLRRRLKSDVALEVPPKREVVVYAPLSKKQEIFYTAIVNRTIAN
MFGSSEKETIELSPTGRPKRRTRKSINYSKIDDFPNELEKLISIQPEVDRERAVEVNIPVESEVNLKL
QNIMMLLRKCCNHPYLIEYPIDPVTQEFKIDEELVTNSGKFLILDRMLPELKKRGHKVLLFSQMTSMLDI
LMDYCHLRDFNFSRLDGSMSYSEREKNMHSFNTDPEVFIFLVSTRAGGLGINLTAADTVIIYDSWNPQS
DLQAQDRCHRIGQTKPVVVYRLVTANTIDQKIVERAAAKRKLKLIHKNHFKGGQSGNLNLSKNFLDPKE
LMELLKSRDYEREIKGSREKVISDKDLELLLDLDRSDLIDQMNASGPIKEKMGIFKILENSEDSSPECLF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	96.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060533
Locus ID:	3070
UniProt ID:	Q9NRZ9
RefSeq Size:	3165
Cytogenetics:	10q23.33
RefSeq ORF:	2514
Synonyms:	ICF4; LSH; Nbla10143; PASG; SMARCA6
Summary:	This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jan 2014]

Product images:



Coomassie blue staining of purified HELLs protein (Cat# [TP312231]). The protein was produced from HEK293T cells transfected with HELLs cDNA clone (Cat# [RC212231]) using MegaTran 2.0 (Cat# [TT210002]).