

Product datasheet for **TP312352L**

Myoglobin (MB) (NM_005368) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human myoglobin (MB), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212352 representing NM_005368 Red =Cloning site Green =Tags(s)
	MGLSDGEWQLVLNVWGKVEADIPGHGQEVLRIRLFGKHPETLEKFDKFKHLKSEDEMKAEDLKKHGATVLTALGGILKKKGHHEAEIKPLAQSHATKHKIPVKYLEFISECIIQVLQSKHPGDFGADAQGAMNKALELFRKDMASNYKELGFQG
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005359
Locus ID:	4151
UniProt ID:	P02144 , A0A1K0FU49
RefSeq Size:	1078



[View online »](#)


Cytogenetics: 22q12.3

RefSeq ORF: 462

Synonyms: PVALB

Summary: This gene encodes a member of the globin superfamily and is predominantly expressed in skeletal and cardiac muscles. The encoded protein forms a monomeric globular haemoprotein that is primarily responsible for the storage and facilitated transfer of oxygen from the cell membrane to the mitochondria. This protein also plays a role in regulating physiological levels of nitric oxide. Multiple transcript variants encoding distinct isoforms exist for this gene. [provided by RefSeq, May 2020]

Product images:



Coomassie blue staining of purified MB protein (Cat# [TP312352]). The protein was produced from HEK293T cells transfected with MB cDNA clone (Cat# [RC212352]) using MegaTran 2.0 (Cat# [TT210002]).