

Product datasheet for **TP312401M**

CDRT15 (NM_001007530) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CMT1A duplicated region transcript 15 (CDRT15), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212401 representing NM_001007530 Red =Cloning site Green =Tags(s) MFSCCFPTSRGCCFRNGGSESLFRRRCRRRLIPHRRRLSPVIRRIQVPQDSLQALAGQATPEIPLGLQL HTVLVQEIQELIEAQTLAGPCAQEVRLPAPAAEPEPAWEEAPPERALELEGAPAKDQTNEELPEITEVP ESIKRRLGRRVPAATPAPRGNLLLQAWMRVHSWASRLFAPNVLPGTGP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001007531</u>
Locus ID:	146822
UniProt ID:	<u>Q96T59</u>
RefSeq Size:	567




[View online »](#)

Cytogenetics: 17p12

RefSeq ORF: 564

Product images:



Coomassie blue staining of purified CDRT15 protein (Cat# [TP312401]). The protein was produced from HEK293T cells transfected with CDRT15 cDNA clone (Cat# [RC212401]) using MegaTran 2.0 (Cat# [TT210002]).