

Product datasheet for **TP312421L**

BTBD17 (NM_001080466) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human BTB (POZ) domain containing 17 (BTBD17), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212421 representing NM_001080466 Red =Cloning site Green =Tags(s)
	MPRRGYSKPGSWGSFWAMLTLVGLVTHAAQRADVGGEEAAGTSINHSQAVLQRLQELLRQGNASDVVLR VQ AAGTDEVRFVFAHRLLLGLHSELFLELLSNQSEAVLQEPQDCAAVFDKFIRYLYCGELTVLLTQAIPLHR LATKYGVSSLQRGVADYMRAGLAGGAGPAVGWYHYAVGTGDEALRESCLQFLAWNLSAVAASTEWSGAVS P ELLWQLLQRSIDLVLQDELELFHALEAWLGRARPPPAVAERALRAIRYPMIPPAQLFQLQARSAALARHGP AVADLLLQAYQFHAASPLHYAKFFDVNGSAFLPRNYLAPAWGAPWWINNPARDRSTSFQTQLGSPGHD A GRRVTWNVLFSPRWLPVSLRPVYADAAGTALPAARPEDGRPRLVTPASSGGDAAGVSFQKTVLVGARQ Q GRLLVRHAYSFHQSSEEAGDFLAHADLQRRNSEYLVENALHLHLIVKPVYHTLIRTPK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	52.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001073935</u>
Locus ID:	388419
UniProt ID:	<u>A6NE02</u>
RefSeq Size:	1678
Cytogenetics:	17q25.1
RefSeq ORF:	1434
Synonyms:	BTBD17A; LGALS3BPL; TANGO10A

Product images:



Coomassie blue staining of purified BTBD17 protein (Cat# [TP312421]). The protein was produced from HEK293T cells transfected with BTBD17 cDNA clone (Cat# [RC212421]) using MegaTran 2.0 (Cat# [TT210002]).