

## Product datasheet for **TP312421M**

### **BTBD17 (NM\_001080466) Human Recombinant Protein**

#### Product data:


|                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|---------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | Recombinant Proteins                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Description:                          | Recombinant protein of human BTB (POZ) domain containing 17 (BTBD17), 100 µg                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Species:                              | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Expression Host:                      | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Expression cDNA Clone or AA Sequence: | >RC212421 representing NM_001080466<br><b>Red</b> =Cloning site <b>Green</b> =Tags(s)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                       | <p>MPRRGYSKPGSWGSEFWAMLTIVGLVTHAAQRADVGGEEAAGTSINHSQAVLQRLQELLRQGNASDVVLR<br/>VQ<br/>AAGTDEVRFVFAHRLLLGLHSELFLELLSNQSEAVLQEPQDCAAVFDKFIRYLYCGELTVLLTQAIPLHR<br/>LATKYGVSSLQRGVADYMRAGLAGGAGPAVGWYHYAVGTGDEALRESCLQFLAWNLSAVAASTEWSGAVS<br/>P<br/>ELLWQLLQRSIDLVLQDELELFHALEAWLGRARPPPAVAERALRAIRYPMIPPAQLFQLQARSAALARHGP<br/>AVADLLLQAYQFHAASPLHYAKFFDVNGSAFLPRNYLAPAWGAPWVINNPARDRSTSFQTLGPGSGHD<br/>A<br/>GRRVTWNVLFSPRWLPVSLRPVYADAAGTALPAARPEDGRPRLVTPASSGGDAAGVSFQKTVLVGARQ<br/>Q<br/>GRLLVRHAYSFHQSSEEAGDFLAHADLQRRNSEYLVENALHLHLIVKPVYHTLIRTPK</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p> |
| Tag:                                  | C-Myc/DDK                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Predicted MW:                         | 52.3 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Concentration:                        | >0.05 µg/µL as determined by microplate BCA method                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Purity:                               | > 80% as determined by SDS-PAGE and Coomassie blue staining                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Buffer:                               | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Preparation:                          | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Note:                                 | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Storage:                              | Store at -80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |



[View online »](#)

|                      |                                                                                                                                               |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Stability:</b>    | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| <b>RefSeq:</b>       | <u>NP_001073935</u>                                                                                                                           |
| <b>Locus ID:</b>     | 388419                                                                                                                                        |
| <b>UniProt ID:</b>   | <u>A6NE02</u>                                                                                                                                 |
| <b>RefSeq Size:</b>  | 1678                                                                                                                                          |
| <b>Cytogenetics:</b> | 17q25.1                                                                                                                                       |
| <b>RefSeq ORF:</b>   | 1434                                                                                                                                          |
| <b>Synonyms:</b>     | BTBD17A; LGALS3BPL; TANGO10A                                                                                                                  |

### Product images:



Coomassie blue staining of purified BTBD17 protein (Cat# [TP312421]). The protein was produced from HEK293T cells transfected with BTBD17 cDNA clone (Cat# [RC212421]) using MegaTran 2.0 (Cat# [TT210002]).