

Product datasheet for **TP312450M**

C20orf151 (RBBP8NL) (NM_080833) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 20 open reading frame 151 (C20orf151), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212450 representing NM_080833 Red =Cloning site Green =Tags(s)

MESFMESLNRLKEIHEKEVLGLQNKLLELNSERCRDAQRIEELFSKNHQLREQQKTLKENLRVLENRLRA
GLCDRCMTQELARKRQQEFESSHLQNLQRIFILTNE MNGLKEENETLKEEVKRLRGLGDRPKPRAKEGT
SDPPSPLLLPSPGGWKAITEKPPGGHEEAEEDHQGVGLRGEEKPAGHRTSPVAKISPGATLPESRAPDMS
PQRISNQLHGTTIAVVRPGSQACPADRG PANGTPPPLPARSSPPSPAYERGLSLDSFLRASRPSAMTHEAP
KLSPKVDRLCLLNRLSLHLQSPHSSPLAPAAAPSDPRLQDLKAREAEAWEEPTELLGLPSALAGMQDLR
LEGALHLLLAQQQLRARARAGSVRPRGQPTPGEMPLPSLPVGS DSEGPENEGTRAALAAAGLSGGRHTQP
A
GPGRAQRTEAAATQDCALDKPLDLSEWGRARGQDTPKPAQHGSLS PAAAHTASPEPPTQSGPLTRSPQ
A
LSNGTKGTRVPEQEEASTPMDPSRPLPGSQLSLSSPGSTEDTGRPLPPHPQPPHPQPPDLGDHPEP
SKAEVLRPEDELDETDTPGSEVGLSSQAEATTSTTGEGPECTQEHGQGPPRKRKRASEPGDKASKKP
SRGRRKLTATEGPGSPRDAEDHSPSPNSSPWEET

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	71.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_543023</u>
Locus ID:	140893
UniProt ID:	<u>Q8NC74</u>
RefSeq Size:	2793
Cytogenetics:	20q13.33
RefSeq ORF:	1992
Synonyms:	C20orf151

Product images:



Coomassie blue staining of purified RBBP8NL protein (Cat# [TP312450]). The protein was produced from HEK293T cells transfected with RBBP8NL cDNA clone (Cat# [RC212450]) using MegaTran 2.0 (Cat# [TT210002]).