

Product datasheet for **TP312451L**

PID1 (NM_017933) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphotyrosine interaction domain containing 1 (PID1), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212451 representing NM_017933 Red =Cloning site Green =Tags(s)
	 MPRIAGNHLMLEESRTCSSPELLDGVWPCQLPHFGLPASEMHFQTMLKSKLNVLTLLKKEPLPAVIFHEPE AIELCTTTPLMKTRTHSGCKVTYLGKVSTTGMQFLSGCTEKPVIELWKKHTLAREDFVPANALLEIRPFQ VWLHHLDDHKGEATVHMDTFQVARIAYCTADHNVSPNIFAWVYREINDDLISYQMDCHAVECESKLEAKKLA HAMMEAFRKTFFHSMKSDGRIHSNSSSSEEVSQELESDDG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	27.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_060403</u>
Locus ID:	55022



[View online »](#)

UniProt ID:	Q7Z2X4
RefSeq Size:	2745
Cytogenetics:	2q36.3
RefSeq ORF:	744
Synonyms:	HMFN2073; NYGGF4; P-CLI1; PCL11
Summary:	Increases proliferation of preadipocytes without affecting adipocytic differentiation. [UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified PID1 protein (Cat# [TP312451]). The protein was produced from HEK293T cells transfected with PID1 cDNA clone (Cat# [RC212451]) using MegaTran 2.0 (Cat# [TT210002]).