

Product datasheet for **TP312596**

SNAP25 (NM_130811) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human synaptosomal-associated protein, 25kDa (SNAP25), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212596 representing NM_130811 Red =Cloning site Green =Tags(s)

MAEDADMRNELEEMQRRADQLADESLESTRMLQLVEESKDAGIRTLVMLDEQGEQLERIEEGMDQINKD
MKEAEKNLTDLGKFCGLCVCPCNKLKSSDAYKKAWGNNQDGVVASQPARVVDEREQMAISGGFIRRVNTD
ARENEMDENLEQVSGIIGNLRHMALDMGNEIDTQNRQIDRIMEKADSNKTRIDEANQRATKMLGSG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	23.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_570824
Locus ID:	6616
UniProt ID:	P60880



[View online »](#)

RefSeq Size: 2053

Cytogenetics: 20p12.2

RefSeq ORF: 618


Synonyms: bA416N4.2; CMS18; dj1068F16.2; RIC-4; RIC4; SEC9; SNAP; SNAP-25; SUP

Summary: Synaptic vesicle membrane docking and fusion is mediated by SNAREs (soluble N-ethylmaleimide-sensitive factor attachment protein receptors) located on the vesicle membrane (v-SNAREs) and the target membrane (t-SNAREs). The assembled v-SNARE/t-SNARE complex consists of a bundle of four helices, one of which is supplied by v-SNARE and the other three by t-SNARE. For t-SNAREs on the plasma membrane, the protein syntaxin supplies one helix and the protein encoded by this gene contributes the other two. Therefore, this gene product is a presynaptic plasma membrane protein involved in the regulation of neurotransmitter release. Two alternative transcript variants encoding different protein isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: SNARE interactions in vesicular transport

Product images:



Coomassie blue staining of purified SNAP25 protein (Cat# TP312596). The protein was produced from HEK293T cells transfected with SNAP25 cDNA clone (Cat# [RC212596]) using MegaTran 2.0 (Cat# [TT210002]).