

Product datasheet for **TP312607M**

MST4 (STK26) (NM_016542) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human serine/threonine protein kinase MST4 (MST4), transcript variant 1, 100 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC212607 protein sequence
Red=Cloning site **Green**=Tags(s)

MAHSPVAVQVPGMQNNIADPEELFTKLERIGKGSFGGEVFKGIDNRTQQVVAIKIIDLEEADEIEDIQQE
 ITVLSQCDSSYVTKYYGKSYLKGSKLWIIMEYLGGSALDLLRAGPFDEFQIATMLKEILKGLDYLHSEKK
 IHRDIKAANVLLSEQGDVKLADFGVAGQLTDTQIKRNTFVGTPFWMAPEVIQQSAYDSKADIWSLGITAI
 ELAKGEPPNSDMHPMRVLFIPKNNPPTLVGDFTESFKEFIDACLNKDPSPRPTAKELLKHKFIVKNSKK
 TSYLTELIDRFKRWKAEGHSDDESSEGSSESTSRNNTHPEWSFTTVRKKPDPKKVQNGAEQDLVQTL
 SCLSMIITPAFAELKQQDENNASRNQAIEELEKSIVAEAAACPGITDKMVKKLIEKFQKCSADESP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 46.3 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_057626](#)




[View online »](#)

Locus ID: 51765
UniProt ID: [Q9P289](#)
RefSeq Size: 3352
Cytogenetics: Xq26.2
RefSeq ORF: 1248
Synonyms: MASK; MST4

Summary: The product of this gene is a member of the GCK group III family of kinases, which are a subset of the Ste20-like kinases. The encoded protein contains an amino-terminal kinase domain, and a carboxy-terminal regulatory domain that mediates homodimerization. The protein kinase localizes to the Golgi apparatus and is specifically activated by binding to the Golgi matrix protein GM130. It is also cleaved by caspase-3 in vitro, and may function in the apoptotic pathway. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Protein Kinase

Product images:



Coomassie blue staining of purified STK26 protein (Cat# [TP312607]). The protein was produced from HEK293T cells transfected with STK26 cDNA clone (Cat# [RC212607]) using MegaTran 2.0 (Cat# [TT210002]).