

Product datasheet for TP312662L

HMBS (NM_001024382) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hydroxymethylbilane synthase (HMBS), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212662 representing NM_001024382 Red=Cloning site Green=Tags(s)

MRVIRVGTRKSQLARIQTDSWATLKASYPGLQFEIIMSTTGDKILDALTALSKIGEKSLFTKELEHALEK
NEVDLVVHSLKDLPTVLPFGFTIGAICKRENPHDAVVFHPKFVGKTLETLPKESVGTSSLRRAAQLQRK
FPHLEFRSIRGNLNLRLRKLDEQQEFSAILATAGLQRMGWVHNRVQILHPEECMYAVGQGALGVEVRAK
DQDILDVLGVLDHPETLLRCIAERAFLRHLEGGCSVPVAVHTAMKDGQLYLTGGVWSLDGSDSIQETMQA
TIHVPAQHEDGPEDDPQLVGITARNIPRGPQLAAQNLGISLANLLSKGAKNILDVARQLNDAH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	37.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001019553
Locus ID:	3145



[View online »](#)


UniProt ID: [P08397](#)
RefSeq Size: 1428
Cytogenetics: 11q23.3
RefSeq ORF: 1032
Synonyms: PBG-D; PBGD; PORC; UPS

Summary: This gene encodes a member of the hydroxymethylbilane synthase superfamily. The encoded protein is the third enzyme of the heme biosynthetic pathway and catalyzes the head to tail condensation of four porphobilinogen molecules into the linear hydroxymethylbilane. Mutations in this gene are associated with the autosomal dominant disease acute intermittent porphyria. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Porphyrin and chlorophyll metabolism

Product images:



Coomassie blue staining of purified HMBS protein (Cat# [TP312662]). The protein was produced from HEK293T cells transfected with HMBS cDNA clone (Cat# [RC212662]) using MegaTran 2.0 (Cat# [TT210002]).