

Product datasheet for **TP312698M**

Histidine decarboxylase (HDC) (NM_002112) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human histidine decarboxylase (HDC), 100 µg

Species: Human

Expression Host: HEK293T

Expression cDNA >RC212698 protein sequence

Clone or AA **Red**=Cloning site **Green**=Tags(s)

Sequence:

MMEPEEYRERGREMVDYICQYLSTVRERRVMPDVQPGYLRAQLPESAPEDPDSWDSIFGDIERIIMPGVV
HWQSPHMHAYYPALTSWPSLLGDMLADAINCLGFTWASSPACTELEMNVMDWLAKMLGLPEHFLHHHPSS
QGGGVLQSTVSESTLIALLAARKNKILEMKTSEPDADDESCLNARLVAYASDQAHSSVEKAGLISLVKMKF
LPVDDNFSLRGEALQKAIEEDKQRGLVPVFCATLGTGVCADFCLSELGPICAREGLWLHIDAAYAGTA
FLCPEFRGFLKGIYADSFTFNPSKWMMVHFDCSTGFWVKDKYKLLQQTFSVNIYLRHANSVATDFMHWQ
IPLSRRFRSVKLVWFVIRSFVGNLQAHVRHGTEMAKYFESLVRNDPSFEIPAKRHLGLVVFRLKGNCLT
ENVLKEIAKAGRLFLIPATIQDKLIIRFTVTSQFTTRDDILRDWNLIRDAATLILSQHCTSQPSPRVGNL
ISQIRGARAWACGTSLSVSGAGDDPVQARKIKQPQRVGVGAPMKRENGLHLETLDPVDDCFSEEAPDA
TKHKLSSFLFSYLSVQTKKKTVRSLSNSVPVSAQKPLPTEASVKNKGSSRVIRFSRFPEDMMMLLKSAF
KKLIKFSVSPSFPECSSQCGLQLPCCPLQAMV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 74 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.


Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.



[View online »](#)

| | |
|--------------------------|---|
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_002103 |
| Locus ID: | 3067 |
| UniProt ID: | P19113 |
| RefSeq Size: | 2646 |
| Cytogenetics: | 15q21.2 |
| RefSeq ORF: | 1986 |
| Summary: | This gene encodes a member of the group II decarboxylase family and forms a homodimer that converts L-histidine to histamine in a pyridoxal phosphate dependent manner. Histamine regulates several physiologic processes, including neurotransmission, gastric acid secretion, inflammation, and smooth muscle tone.[provided by RefSeq, Aug 2010] |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Histidine metabolism, Metabolic pathways |

Product images:

Coomassie blue staining of purified HDC protein (Cat# [TP312698]). The protein was produced from HEK293T cells transfected with HDC cDNA clone (Cat# [RC212698]) using MegaTran 2.0 (Cat# [TT210002]).