

Product datasheet for **TP312782M**

KLHL17 (NM_198317) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human kelch-like 17 (Drosophila) (KLHL17), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212782 representing NM_198317 Red =Cloning site Green =Tags(s)

MQPRSERPAGRTQSPEHGSPGPGPEAPPPPPQPPAPEAERTRPRQARPAAPMEGAVQLLSREGHSAH
N
SKRHYHDAFVAMSRMRQRGLLCDIVLHVAAKEIRAHKVVLASCSFYHAMFTNEMSESQRQTHVTLHDIDP
QALDQLVQFAYTAEIVVGEGNVQTLLPAASLLQLNGVRDACCKFLLSQLDPSNCLGIRGFADAHSCDLL
KAAHRYVLQHFDVAKTEEFMLLPLKQVLELVSSDSLNPSEEEVYRAVLSWVKHDVDARRQHVPRLMKC
VRLPLLSRDFLLGHVDAESLVRHHPDCKDLLIEALKFHLLPEQRGVLGTSRTRPRRCEGAGPVLFAVGGG
SLFAIHGDCEAYDTRTDRWHVASMSTRRARVGVAAVGNRLYAVGGYDGTSDLATVESYDPVTNTWQPE
V
SMGTRRSCLGVAALHGLLYSAGGYDGASCLNSAERYDPLTGTWTSVAAMSTRRRYVRVATLDGNLYAVGG
YDSSSHLATVEKYEPQNVWSPVASMLSRSSAGVAVLEGALYVAGGNDGTSCLSNVERYSPKAGAWESV
APMNIRRSTHDLVAMDGWLYAVGGNDGSSSLNSIEKYNPRTNKWVAASCMFTRRSSVGAVLELLNFPP
P
SSPTLSVSSTSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	69.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_938073
Locus ID:	339451
UniProt ID:	Q6TDP4
RefSeq Size:	2560
Cytogenetics:	1p36.33
RefSeq ORF:	1926
Synonyms:	AF
Summary:	The protein encoded by this gene is expressed in neurons of most regions of the brain. It contains an N-terminal BTB domain, which mediates dimerization of the protein, and a C-terminal Kelch domain, which mediates binding to F-actin. This protein may play a key role in the regulation of actin-based neuronal function. [provided by RefSeq, Aug 2010]

Product images:



Coomassie blue staining of purified KLHL17 protein (Cat# [TP312782]). The protein was produced from HEK293T cells transfected with KLHL17 cDNA clone (Cat# [RC212782]) using MegaTran 2.0 (Cat# [TT210002]).