

Product datasheet for TP312795L

ACBD5 (NM_145698) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human acyl-Coenzyme A binding domain containing 5 (ACBD5), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212795 representing NM_145698 Red =Cloning site Green =Tags(s)
	MADTRSVHETRFEAAVKVIQSLPKNGSFQPTNEMMLKFYSFYKQATEGPCKLSRPGFWDPIGRYKWDAWS SLGDMTKEEAMIAYVEEMKKIETMPMTEKVEELLRIGPFYEIVEDKKSGRSSDITSDLGNVLTSTPNA KTVNGKAESSDGAEESEEEEAQEEVKGAEQSDNDKMMKKSADHKNLEVITNGYDKDGFVQDIQNDIHA SSSLNGRSTEEVKPIDENLGQTGKSAVCIHQDINDDHVEDVTGIQHLSDSSEVYCDSEMEQFGQEEESLD SFTSNNGPFQYYLGGHSSQPMENSGFREDIQVPPGNGNIGNMQVWAVEGKGEVKHGGEDGRNNSGAPHRE KRGGETDEFSNRRRGRGRMQHLSEGTKGRQVGGSDGERWGS DRGSRGSLNEQIALVLMRLQEDMQNVL QRLQKLETLTALQAKSSTSTLQTAPQPTSQRPSWWPFEMSPGVLTFAIWPFAIQWLVLVYLYQRRRRKLN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	54.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_663736](#)

Locus ID: 91452

UniProt ID: [Q5T8D3](#), [B7Z2A7](#), [Q8NCM9](#)

RefSeq Size: 2327

Cytogenetics: 10p12.1


RefSeq ORF: 1470

Synonyms: RDLKD

Summary: This gene encodes a member of the acyl-Coenzyme A binding protein family, known to function in the transport and distribution of long chain acyl-Coenzyme A in cells. This gene may play a role in the differentiation of megakaryocytes and formation of platelets. A related protein in yeast is involved in autophagy of peroxisomes. A mutation in this gene has been associated with autosomal dominant thrombocytopenia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified ACBD5 protein (Cat# [TP312795]). The protein was produced from HEK293T cells transfected with ACBD5 cDNA clone (Cat# [RC212795]) using MegaTran 2.0 (Cat# [TT210002]).