

Product datasheet for TP312799M

NCRNA00174 (LINC00174) (NM_181722) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical protein LOC285908 (LOC285908), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212799 representing NM_181722 Red =Cloning site Green =Tags(s) MASAGPNRPEISLARNSTCVGCPNNHSFSRTVAHGGALASAWPPQAQFLPVGGRSRPGSCPSASSPGP E LVPVGLSRPSSPGHLHGPSSCLTTTTFGPAPAQLLAAVVGPRLPCVQASRTHLRISGGPERPGSCLPAAS PGAAASRRPPQATFPASRQPRQARLPPAGGLRLRLISCPAAASPGQAPACRQAPQAQLLRPEGFSRPG SCLAAAS TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	21.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_859073</u>
Locus ID:	285908


[View online »](#)

RefSeq Size: 4429
 Cytogenetics: 7q11.21
 RefSeq ORF: 651
 Synonyms: NCRNA00174

Product images:



Coomassie blue staining of purified LOC285908 protein (Cat# [TP312799]). The protein was produced from HEK293T cells transfected with LOC285908 cDNA clone (Cat# [RC212799]) using MegaTran 2.0 (Cat# [TT210002]).