

Product datasheet for **TP312833L**

TOGARAM2 (NM_199280) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 179, member A (FAM179A), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212833 representing NM_199280 Red =Cloning site Green =Tags(s)

MGTRDDVPEAKVLVPVAVYCGSIPRTSAGPRVLPPGSINSSLPHGEGSLQPEPRALLNNEEPSQLLRGLG
QLGGLKLDTPSKGWQARNGHPRNLRALS LGDQPLVLLSPSEANSVARDTIQIKDKLKKRRLSEGLAAS
SRASLDPGGGPQGVPLHSTIPRATSQRLLRVPRMPLIQSIPTTPEASGVKEKGLDLP GSIPGPHEL RP G
AQEAQISWQYLHCNDEKMQKSLGAIVIPPIPKARTVAATPSRVPGSLPSPLPPGQGVLTGLRAPRTRLAR
GSGPREKTPASLEPKPLASPIRDRPAAAKKPALPFSQSAPTTLTAFSFD CAREACPLKEEDQKEIGTKIQ
VTISKSAREKMQLKQMKEMELLRRLEERTGQELTSQCLGSQRAFMKEGLPLRGSGT LSVPTRLSGPCR
NDVSIILRKWASRASLPSIPISRQEPFRFARHASANSLPAVLT LGSPWEEEEEMDLRACKELRPFSNP EL
GLRDALQCLNSSDWQMKEKGLVSIQRLAACHSEVLTGKLHDVCLVVTGEVTNLR SKVSHLAISTLGLDFQ
ALKKNMDQEAEEIARCLLQK MADTNEFIQRAAGQSLRAMVENVT LARSLVLT SAGVYHRNPLIRKYAAE
HLSAVLEQIGA EKLLSGTRDSTDMLVHNLVRLAQDSNQDTRFYGRKMVNILMANTKFDAFLKQSLPSYDL
QKVMAAIKQQGIEDNDELPSAKGRKVLRLVWCENGLPIKEGLSCNGPRLVGLRSTLQGRGEMVEQLREL
TRLLEAKDFRSRM EGVGQ LLELCKAKTELVT AHLVQVFDAFTPRLQDSNKKVNQWALESF AKMIPLLRES
LHPMLLSIIITVADNLNSKN SGIYAAAVLDAMVESLDNLCLLPALAGRVRFLSGRAVL DVTDRLAVLV
ASVYPRKPQ AVERHVLPI LWHFLNTATRNGTLPGPSGNIRGVCRLSRSLQE HMGSRLLDFAASQPKHVL
KTLQELLDSESLGGSRKATDRGVAPDSKTTGSSYPFQLD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	111 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_954974
Locus ID:	165186
UniProt ID:	Q6ZUX3
RefSeq Size:	3564
Cytogenetics:	2p23.2
RefSeq ORF:	3057
Synonyms:	FAM179A

Product images:



Coomassie blue staining of purified TOGARAM2 protein (Cat# [TP312833]). The protein was produced from HEK293T cells transfected with TOGARAM2 cDNA clone (Cat# [RC212833]) using MegaTran 2.0 (Cat# [TT210002]).