

## Product datasheet for **TP312850L**

### hnRNP F (HNRNPF) (NM\_001098204) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens heterogeneous nuclear ribonucleoprotein F (HNRNPF), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212850 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MMLGPEGGEGFVVKLRGLPWSCSVEDVQNFLSDCTIHDGAAGVHFIYTREGRQSGEAFVELGSEDDVKMALKKDRESMGHRYIEVFRSHRTEMDWVLKHSGPNSADSANDGFVRLRGLPFGCTKEEIVQFFSGLIEVPNGITLPVDPEGKITGEAFVQFASQELAELKALGKHKERIGHRYIEVFKSSQEEVRSYSDPPLKFMVQRP GPYDRPGTARRYIGIVKQAGLERMRPGAYSTGYGGYEEYSGLSGDYGFITDLFGRDLSYCLSGMYDHR YGDSEFTVQSTTGHCVHMRGLPYKATENDIYNFFSPLNPVRVHIEIGPDGRVTGEADVEFATHEEAVAAMSKDRANMQHRYIELFLNSTTGASNGAYSSQVMQGMGVSAQAATYSGLESQSVSGCYGAGYSGQNSMGGYD

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	45.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001091674</a>




[View online »](#)

Locus ID:	3185
UniProt ID:	<a href="#">P52597</a> , <a href="#">A0A024R7T3</a>
RefSeq Size:	2650
Cytogenetics:	10q11.21
RefSeq ORF:	1245
Synonyms:	HNRPF; mcs94-1; OK/SW-cl.23

**Summary:** This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and regulate alternative splicing, polyadenylation, and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNAs which have guanosine-rich sequences. This protein is very similar to the family member hnRPH. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified HNRNPF protein (Cat# [TP312850]). The protein was produced from HEK293T cells transfected with HNRNPF cDNA clone (Cat# [RC212850]) using MegaTran 2.0 (Cat# [TT210002]).