

Product datasheet for **TP312899L**

AGXT (NM_000030) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human alanine-glyoxylate aminotransferase (AGXT), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC212899 protein sequence Red =Cloning site Green =Tags(s)

MASHKLLVTPPKALLKPLSIPNQLLLGPGPSNLPPRIMAAGGLQMIGSMSKDMYQIMDEIKEGIQYVFQT
RNPLTLVISGSGHCALEAALVNVLEPGDSFLVGANGIWGQRAVDIGERIGARVHPMTKDPGGHYTLQEVE
EGLAQHKPVLLFLTHGESSTGVLQPLDGFELCHRYKCLLLVDSVASLGGTPLYMDRQIDILYSGSQKA
LNAPPGTSLISFSDKAKKKMYSRKT KPFSFYLDIKWLANFWGCDDQPRMYHHTIPVISLYSLRESLALIA
EQGLENSWRQHREAAAYLHGRLQALGLQLFVKDPALRLPTVTTVAVPAGYDWRDIVSYVIDHFDIEIMGG
LGPSTGKVLRIPLLGCNATRENVDRVTEALRAALQHCPKKKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	42.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_000021</u>
Locus ID:	189



[View online »](#)

UniProt ID:	P21549
RefSeq Size:	1611
Cytogenetics:	2q37.3
RefSeq ORF:	1176
Synonyms:	AGT; AGT1; AGXT1; PH1; SPAT; SPT; TLH6
Summary:	This gene is expressed only in the liver and the encoded protein is localized mostly in the peroxisomes, where it is involved in glyoxylate detoxification. Mutations in this gene, some of which alter subcellular targeting, have been associated with type I primary hyperoxaluria. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Alanine, aspartate and glutamate metabolism, Glycine, serine and threonine metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified AGXT protein (Cat# [TP312899]). The protein was produced from HEK293T cells transfected with AGXT cDNA clone (Cat# [RC212899]) using MegaTran 2.0 (Cat# [TT210002]).