

Product datasheet for TP313082

Transglutaminase 7 (TGM7) (NM_052955) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human transglutaminase 7 (TGM7), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213082 representing NM_052955 Red=Cloning site Green=Tags(s)
	<p>MDQVATLRLESVDLQSSRNKEHHTQEMGVKRLTVRRGQPFYLRLSFSRPFQSQNDHITFVAETGPKPSE LLGTRATFFLTRVQPGNVWSASDFTIDSNSLQVSLFTPANAVIGHYTLKIEISQGGHSTYPLGTFILL FNPWSPEDDVYLPSEILLQEYIMRDYGFVYKGERFITSWPWNYGQFEEDIIDICFEILNKSPLYHLKNPA KDCSQRNDVVYVCRVVSAMINSNDDNGVLQGNWGEDYSKGVSPLEWKGSVAIQQWSARGGQPVKYGQCW VFASVMCTVMRCLGVPTRVVSNFRSAHNVDNRNLTIIDYDRNAEMLSTQKRDKIWNFHVWNECWMIRKDL PPGYNGWQVLDPTPQQTSSGLFCCGPASVKAIREGDVHLAYDTPFVYAENVNADEVIWLLGDGQAQEILAH NTSSIGKEISTKMVGSDQRQSITSSYKYPEGSPEERAVFMKASRKMLGPPRSLPFLDLLESGGLRDQPA QLQLHLARPEWGQDLQLLLRIQRPDSTHPRGPIGLVVRFCAQALLHGGGTQKPFWRHTVRMNLDFGKE TQWPLLLPYSNYRNKLTDEKLIRVSGIAEVEETGRSMLVLKIDICLEPPHLSIEVSERAIEVKGALRVHVTL TNTLMVALSSCTMVLEGSGLINGQIAKDLGTLVAGHTLQIQLDLYPTKAGPRQLQVLISSNEVKEIKGYK DIFVTVAGAP</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	79.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_443187](#)

Locus ID: 116179

UniProt ID: [Q96PF1](#)

RefSeq Size: 2312

Cytogenetics: 15q15.2-q15.3


RefSeq ORF: 2130

Synonyms: TGMZ

Summary: Transglutaminases (TGM; EC 2.3.2.13) are a family of structurally and functionally related enzymes that stabilize protein assemblies through the formation of gamma-glutamyl-epsilon lysine crosslinks. For additional background information on transglutaminases, see TGM1 (MIM 190195). [supplied by OMIM, Jul 2002]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified TGM7 protein (Cat# TP313082). The protein was produced from HEK293T cells transfected with TGM7 cDNA clone (Cat# [RC213082]) using MegaTran 2.0 (Cat# [TT210002]).