

Product datasheet for **TP313133M**

Midkine (MDK) (NM_002391) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human midkine (neurite growth-promoting factor 2) (MDK), transcript variant 3, 100 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC213133 representing NM_002391
Red=Cloning site **Green**=Tags(s)

MQHRGFLLLTLLALLALTSAVAKKKDKVKKGGPGSECAEWAWGPCTPSSKDCGVGFREGTCGAQTQRIRC
RVPCNWKKEFGADCKYKFENWGACDGGTGTKVRQGLKKARYNAQCQETIRVTKPCTPKTKAKAKAKKGGK
GKD

TRTRPLE**Q**KLISEEDLAANDILDY**K**DDDD**KV**

Tag: C-Myc/DDK

Predicted MW: 15.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_002382](#)

Locus ID: 4192

UniProt ID: [P21741](#)



[View online »](#)

RefSeq Size: 866

Cytogenetics: 11p11.2


RefSeq ORF: 429

Synonyms: ARAP; MK; NEGF2

Summary: This gene encodes a member of a small family of secreted growth factors that binds heparin and responds to retinoic acid. The encoded protein promotes cell growth, migration, and angiogenesis, in particular during tumorigenesis. This gene has been targeted as a therapeutic for a variety of different disorders. Alternatively spliced transcript variants encoding multiple isoforms have been observed. [provided by RefSeq, Jul 2012]

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Product images:



Coomassie blue staining of purified MDK protein (Cat# [TP313133]). The protein was produced from HEK293T cells transfected with MDK cDNA clone (Cat# [RC213133]) using MegaTran 2.0 (Cat# [TT210002]).