

Product datasheet for TP313157

MEGF6 (NM_001409) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Homo sapiens multiple EGF-like-domains 6 (MEGF6), 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA >RC213157 representing NM_001409

Clone or AA Red=Cloning site Green=Tags(s)

Sequence:

MSFLEEARAAGRAVVLALVLLLLPAVPVGASVPPRLLPLQPGMPHVCAEQELTLVGRRQPCVQALSHTV
PVWKAGCGWQAWCVGHERRTVYYMGYRQVYTTEARTVLRCRGMWQPDDEEGCLSAECSASLCFHGGRCV
PGSAQPCHCPPGFQGPCQYDVDECRTNNGGCQHRCVNTPGSYLCECKPGFRLHTDSRTCLAINSCALGN
GGCQHHCVQLTITRHCRCRPGFQLQEDGRHCVRRSPCANRNGSCMHRCQVVRGLARCECHVGYQLAADG
KACEDVDECAAGLAQCAHGCLNTQGSFKCVCHAGYELGADGRQCYRIEMEIVNSCEANNGGCSHGCSHTS
AGPLCTCPRGYELDTDQRTCIDVDDCADSPCCQVCTNPNPGGYECGYAGYRLSADGCGCEDVDECASSR
GGCEHHCTNLAGSFQCSCEAGYRLHEDRRGCSPLEEPMVDLDGELPFVRPLPHIAVLQDELPLQFQDDD
GADEEEAELRGEHTLTEKFVCLDDSFHDCSLTCDDCRNGGTCLLGLDGCDCEPWTGLICNETCPPDTF
GKNCSFSCSCQNGGTCDSVTGACRCPPGVSGTNCEDGCPKGYGKHCRKCCNCANRGRCHRLYGACLCDP
GLYGRFCHLTCPWAFPGCSEECQCVQPHQTQSCDKRDGSCSCKAGFRGERCQAECELGYPGPGCWQACT
CPVGACDSVSGECGKRCPAGFQGEDCGQECVPVGTGFGVNCSSSCSGGAPCHGVTGQCRCPPGRTGEDCE
ADCPEGRWGLGCQEICPACQHAARCDPETGACLCLPGFVGSRCQDVCPAGWYGPSCQTRCSCANDGHCHP
ATGHCSCAPGWTGFSCQRACDTGHWGPDCSHPCNCSAGHGSCDAISGLCLCEAGYVGPCEQQCPQGHFG
PGCEQRCQCQHGAACDHVSGACTCPAGWRGTFCHEACPAFFGLDCRSACNCTAGAACDAVNGSCLCPAG
RRGPRCAETCPAHTYGHNCQACACFNGASCDPVHGQCHCAPGWMGPSCQACPAAGLYGDNCRHSCLCQN
GGTCDPVSGHCACPEGWAGLACEKECLPRDVRAGCRHSGGCLNGLCDPHTGRCLCPAGWTGDKCQSPCL
RGWFGCAQRCSCPPGAACHHVTGACRCPPGFTGSGCEQACPPGSFGEDCAQMCQCPGENPACHPATGT
CSCAAGYHGPSCQRCPPGRYGPGEQLCGCLNNGSCDAATGACRCPTGFLGTDNLTCPQGRFGPNCTH
VCGCGQGAACDPVTGTCLCPPGRAGVRCERGCQNRFGVGCEHTCSCRNGGLCHASNGSCSCLGWTGRH
CELACPPGRYGAACHLECSCHNNSTCEPATGTCRCGPGFYGQACEHPCPPGFHAGAGCQGLWCQHGAAPCD
PISGRCLCPAGFHGHFCERGCEPGSFGECHQRCDGCDGAPCDPVTGLCLCPPGRSGATCNLDCRRGQFG
PSCTLHCDCGGGADCDPVSGQCHCVDGYMGPTCREGGPLRLPENPSLAQGSAGTLPASSRPTSRSGGPAR
H


TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Tag:	C-Myc/DDK
Predicted MW:	161 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001400</u>
Locus ID:	1953
UniProt ID:	<u>O75095</u>
RefSeq Size:	7450
Cytogenetics:	1p36.32
RefSeq ORF:	4623
Synonyms:	EGFL3
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified MEGF6 protein (Cat# TP313157). The protein was produced from HEK293T cells transfected with MEGF6 cDNA clone (Cat# [RC213157]) using MegaTran 2.0 (Cat# [TT210002]).