

Product datasheet for **TP313276M**

PRRX1 (NM_006902) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human paired related homeobox 1 (PRRX1), transcript variant pmx-1a, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213276 representing NM_006902 Red =Cloning site Green =Tags(s)

MTSSYGHVLERQPALGGRLDSPGNLDTLQAKKNFSVSHLLDLEEAGDMVAAQADENVGEAGRSLLLESPGL
TSGSDTPQQDNDQLNSEKKKRKQRRNRRTTFNSSLQALERVFERHYPDAFVREDLARRVNLTEARVQV
WFQNRRAKFRNRERAMLANKNASLLKSYSGDVTAVEQPIVPRPAPRPTDYLSWGTASPYRSSSLPRCCLH
EGLHNGF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	24.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_008833</u>
Locus ID:	5396




[View online »](#)

UniProt ID: [P54821](#)
RefSeq Size: 4071
Cytogenetics: 1q24.2
RefSeq ORF: 651
Synonyms: AGOTC; PHOX1; PMX1; PRX-1; PRX1

Summary: The DNA-associated protein encoded by this gene is a member of the paired family of homeobox proteins localized to the nucleus. The protein functions as a transcription co-activator, enhancing the DNA-binding activity of serum response factor, a protein required for the induction of genes by growth and differentiation factors. The protein regulates muscle creatine kinase, indicating a role in the establishment of diverse mesodermal muscle types. Alternative splicing yields two isoforms that differ in abundance and expression patterns. [provided by RefSeq, Jul 2008]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified PRRX1 protein (Cat# [TP313276]). The protein was produced from HEK293T cells transfected with PRRX1 cDNA clone (Cat# [RC213276]) using MegaTran 2.0 (Cat# [TT210002]).