

Product datasheet for TP313375M

Caspase 9 (CASP9) (NM_032996) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human caspase 9, apoptosis-related cysteine peptidase (CASP9), transcript variant beta, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213375 representing NM_032996 Red=Cloning site Green=Tags(s)
	MDEADRRLLRRCLRLVEELQVDQLWDVLLSRELFPHMIEDIQRAGSGSRRDQARQLIIDLETRGSQAL PLFISCLEDTGQDMLASFLRTNRQAAKLSKPTLENLTPVVL RPEIRKPEVLRPETPRPVDIGSGGFGDVE QKDHGFEVASTSPEDESPGSNPEPDATPFQEGLRTFDQLDAISSLPSPDIFVSYSTFPGFVSWRDPKSG SWYVETLDDIFEQWAHSEDLQSLLLRVANAVSVKGIYKQMPGCFNFLRKKLFFKTS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	30 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_127463
Locus ID:	842



[View online »](#)


UniProt ID:	P55211
RefSeq Size:	1584
Cytogenetics:	1p36.21
RefSeq ORF:	798
Synonyms:	APAF-3; APAF3; ICE-LAP6; MCH6; PPP1R56

Summary: This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein can undergo autoproteolytic processing and activation by the apoptosome, a protein complex of cytochrome c and the apoptotic peptidase activating factor 1; this step is thought to be one of the earliest in the caspase activation cascade. This protein is thought to play a central role in apoptosis and to be a tumor suppressor. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013]

Protein Families: Druggable Genome, Protease, Stem cell - Pluripotency

Protein Pathways: Alzheimer's disease, Amyotrophic lateral sclerosis (ALS), Apoptosis, Colorectal cancer, Endometrial cancer, Huntington's disease, Non-small cell lung cancer, p53 signaling pathway, Pancreatic cancer, Parkinson's disease, Pathways in cancer, Prostate cancer, Small cell lung cancer, VEGF signaling pathway, Viral myocarditis

Product images:



Coomassie blue staining of purified CASP9 protein (Cat# [TP313375]). The protein was produced from HEK293T cells transfected with CASP9 cDNA clone (Cat# [RC213375]) using MegaTran 2.0 (Cat# [TT210002]).