

Product datasheet for TP313426

C22orf9 (KIAA0930) (NM_015264) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 22 open reading frame 9 (C22orf9), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213426 representing NM_015264 Red=Cloning site Green=Tags(s)
	MGSQAAAEWRNWASWEVSSSLSGCSMGCFKDDRIVFWTWMFSTYFMEKWAPRQDDMLFYVRRKLAYSGSE SGADGRKAAEPEVEVEVYRRDSKKLPGLGDPDIDWEESVCLNLILQKLDYMTCAVCTRADGGDIHIHKK KSQQVFASPSKHPMDSKGEEKISYPNIFFMIDSFEVFSDMTVGEGEMVCVELVASKTNTFQGVIFQG SIRYEALKKVYDNRVSAARMAQKMSFGFYKYSNMEFVRMKGPGKGHAEMAVSRVSTGDTSPCGTEEDS SPASPMHERVTSFSTPPTPERNNRPAFFPSLKRKVPNRNIAEMKKSHSANDSEEFFREDDGGADLHNAT NLRSRSLSGTGRSLVGSWLKLNLRADGNFLLYAHLTYVTLPLHRILTDILEVRQKPILMT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	46 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_056079</u>



[View online »](#)

Locus ID: 23313
 UniProt ID: [Q6ICG6](#)
 RefSeq Size: 6287
 Cytogenetics: 22q13.31
 RefSeq ORF: 1227
 Synonyms: C22orf9; LSC3

Product images:


Coomassie blue staining of purified KIAA0930 protein (Cat# TP313426). The protein was produced from HEK293T cells transfected with KIAA0930 cDNA clone (Cat# [RC213426]) using MegaTran 2.0 (Cat# [TT210002]).