

Product datasheet for **TP313426L**

C22orf9 (KIAA0930) (NM_015264) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 22 open reading frame 9 (C22orf9), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213426 representing NM_015264 Red =Cloning site Green =Tags(s)
	MGSQAAAEWRNWASWEVSSSLSGCSMGCFKDDRIVFWTWMFSTYFMEKWAPRQDDMLFYVRRKLAYS GSE SGADGRKAAEPEVEVEVYRRDSKKLPGLGDPDIDWEESVCLNLILQKLDYMTCAVCTRADGGDIHHKK KSQQVFASPSKHPMDSKGEESKISYPNIFMIDSFEEVFSDMTVGEEMVCVELVASDKTNTFQGVIFQG SIRYEALKKVYDNRVSVAARMAQKMSFGFYKYSNMEFVRMKGPQGKGHAEMAVSRVSTGDTSPCGTEED S SPASPMHERVTSFSTPPTPERNNRPAFFSPSLKRKVPRNRIAEMKKSHSANDSEEFFREDDGGADLHNAT NLRSRSLSGTGRSLVGSWLKLNRADGNFLLYAHLYVTLP LHRILTDILEVRQKPILMT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	46 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_056079</u>
Locus ID:	23313
UniProt ID:	<u>Q6ICG6</u>
RefSeq Size:	6287
Cytogenetics:	22q13.31
RefSeq ORF:	1227
Synonyms:	C22orf9; LSC3

Product images:



Coomassie blue staining of purified KIAA0930 protein (Cat# [TP313426]). The protein was produced from HEK293T cells transfected with KIAA0930 cDNA clone (Cat# [RC213426]) using MegaTran 2.0 (Cat# [TT210002]).