

Product datasheet for **TP313538M**

AKR1C2 (NM_001354) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human aldo-keto reductase family 1, member C2 (dihydrodiol dehydrogenase 2; bile acid binding protein; 3-alpha hydroxysteroid dehydrogenase, type III) (AKR1C2), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213538 representing NM_001354 Red =Cloning site Green =Tags(s)
	MDSKYQCVKLNLDGHFMPVLFVGFYAPAEVPKSKALEAVKLAIEAGFHHIDSAHVYNNEEQVGLAIRSKIA DGSVKREDIFYTSKLWSNSHRPELVRPALERSLKNLQLDYVDLYLIHFPVSVKPGEEVIPKDENGKILFD TVDLWCATWEAMEKCKDAGLAKSIGVSNFNHRLLEMILNKPGLKYKPVNCQVECHPYFNQRKLLDFCKSKD IVLVAYSALGSHREEPWVDPNSPVLLEDVPLCALAKKHKRTPALIALRYQLQRGVWVLAQSYNEQRIRQN VQVFEFQLTSEEMKAIDGLNRNVRYLTLDFAGPPNYPFSDEY
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	36.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001345



[View online »](#)

Locus ID:	1646
UniProt ID:	P52895
RefSeq Size:	1663
Cytogenetics:	10p15.1
RefSeq ORF:	969
Synonyms:	AKR1C-pseudo; BABP; DD; DD-2; DD/BABP; DD2; DDH2; HAKRD; HBAB; MCDR2; SRXY8; TDD
Summary:	This gene encodes a member of the aldo/keto reductase superfamily, which consists of more than 40 known enzymes and proteins. These enzymes catalyze the conversion of aldehydes and ketones to their corresponding alcohols using NADH and/or NADPH as cofactors. The enzymes display overlapping but distinct substrate specificity. This enzyme binds bile acid with high affinity, and shows minimal 3-alpha-hydroxysteroid dehydrogenase activity. This gene shares high sequence identity with three other gene members and is clustered with those three genes at chromosome 10p15-p14. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]
Protein Families:	Druggable Genome
Protein Pathways:	Metabolism of xenobiotics by cytochrome P450

Product images:



Coomassie blue staining of purified AKR1C2 protein (Cat# [TP313538]). The protein was produced from HEK293T cells transfected with AKR1C2 cDNA clone (Cat# [RC213538]) using MegaTran 2.0 (Cat# [TT210002]).