

## Product datasheet for TP313556M

### TIAM2 (NM\_001010927) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens T-cell lymphoma invasion and metastasis 2 (TIAM2), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213556 representing NM_001010927 Red=Cloning site Green=Tags(s)

MEGPRENQDPPPRPLARHLSADRLRKVIQELVDTEKSYVKDLSCLFELYLEPLQNETFLTQDEMESLFG  
SLPEMLEFQKVFLLETLEDGISASSDFNTLETPSQFRKLLFSLGGSFLLYADHFKLYSGFCANHIKVVQKVL  
ERAKTDKAFKAFLDARNPTKQHSSTLESYLIKPVQRVLKYPLLLKELVSLTDQESEEHYHLTEALKAMEK  
VASHINEMQKIYEDYGTVFDQLVAEQSGTEKEVTELSMGELLMHSTVSWLNPFLSLGKARKDLELTVFVF  
KRAVILVYKENCKLKKKLPNSRPAHNSTDLDPFKFRWLIPISALQVRLGNPAGTENNSIWELIHTKSEI  
EGRPETFQLCCSDSESKTNIVKVIRSIKRENFRRIKCELPLEKTCKDRLVPLKNRVPVSAKLASSRSL  
KVLKNSSSNEWTGETGKGTLLDSDDEGLSSGTQSSGCPTAEGRQDSKSTSPGKYPHPGLADFADNLIKES  
DILSDEDDHRQTVKQGSPTKDIEIQFQRLRISEDPDVHPEAEQPGPESGEGQKGGEQPKLVGRGHFCPI  
KRKANSTKRDRGTLLKAQIRHQSLDSQSENATIDLNSVLEREFVQSLTSVWSEECFYETESHGKS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	70.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_001010927](#)

**Locus ID:** 26230

**UniProt ID:** [Q8IVF5](#)

**RefSeq Size:** 2668

**Cytogenetics:** 6q25.2-q25.3

**RefSeq ORF:** 1878


**Synonyms:** STEF; TIAM-2

**Summary:** This gene encodes a guanine nucleotide exchange factor. A highly similar mouse protein specifically activates ras-related C3 botulinum substrate 1, converting this Rho-like guanosine triphosphatase (GTPase) from a guanosine diphosphate-bound inactive state to a guanosine triphosphate-bound active state. The encoded protein may play a role in neural cell development. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome

**Protein Pathways:** Chemokine signaling pathway, Regulation of actin cytoskeleton

### Product images:



Coomassie blue staining of purified TIAM2 protein (Cat# [TP313556]). The protein was produced from HEK293T cells transfected with TIAM2 cDNA clone (Cat# [RC213556]) using MegaTran 2.0 (Cat# [TT210002]).