

Product datasheet for **TP313671M**

TP53TG3 (NM_016212) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human TP53 target 3 (TP53TG3), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213671 protein sequence Red =Cloning site Green =Tags(s)
	MRASPCISQPAASWHPRPSALRPTAGSGPDTRTPGTVEDGSAPCPAFRSPAVSPCGEPPCCFQISPAEET LELGRLVSPGNCDTLPRAAGFYACHVRSRSLIPCRSTKGRWPLTASAAGLSSFSG
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	12.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057296
Locus ID:	24150
UniProt ID:	Q9ULZ0
RefSeq Size:	1431
Cytogenetics:	16p11.2




[View online »](#)

RefSeq ORF: 372

Synonyms: P53TG3; TP53TG3A; TP53TG3E; TP53TG3F

Product images:



Coomassie blue staining of purified TP53TG3 protein (Cat# [TP313671]). The protein was produced from HEK293T cells transfected with TP53TG3 cDNA clone (Cat# [RC213671]) using MegaTran 2.0 (Cat# [TT210002]).