

Product datasheet for **TP313680**

CCDC90B (NM_021825) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human coiled-coil domain containing 90B (CCDC90B), 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC213680 representing NM_021825

Red=Cloning site **Green**=Tags(s)

MNSRQAWRLFLSQGRGDRWVSRPRGHFSPALRREFFTTTTKEGYDRRPVDITPLEQRKLTFDTHALVQDL
ETHGFDKTQAETIVSALTALSNVSLDTIYKEMVTQAQQEITVQQLMAHLDAIRKDMVILEKSEFANLRAE
NEKMKIELDQVKQQLMHETSRIRADNKLDINLERSRVTD MFTDQEKQLMETTTEFTKKDTQTKSISETS
NKIDAEIASLKTLMESNKLETIRYLAASVFTCLAIALGFYRFWK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 29.3 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_068597](#)


Locus ID: 60492

UniProt ID: [Q9GZT6](#)



[View online »](#)

RefSeq Size:	1585
Cytogenetics:	11q14.1
RefSeq ORF:	762
Synonyms:	MDS011; MDS025
Protein Families:	Transmembrane

Product images:

Coomassie blue staining of purified CCDC90B protein (Cat# TP313680). The protein was produced from HEK293T cells transfected with CCDC90B cDNA clone (Cat# [RC213680]) using MegaTran 2.0 (Cat# [TT210002]).