

Product datasheet for **TP313736**

FAM9B (NM_205849) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human family with sequence similarity 9, member B (FAM9B), 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC213736 representing NM_205849
Red=Cloning site **Green**=Tags(s)

MAAWGKKHAGKDPVRDECEERNRFTETREEDVTDEHGEREPFAETDEHTGANTKKPEDTAEDLTAKRKR
KMDKTCSTKNKSKHALRKKQLKRQKRDYIHSKLLNVLEEYITDEQKEEEEEEGEEEEELIRIFQEQQKK
WQQYRSVRRERLKEMKLLRDQFVKALEDFEDLCDRVFSDEDESELDN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 22.3 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_995321](#)
Locus ID: 171483
UniProt ID: [Q8IZU0](#), [A0A024RBV3](#)
RefSeq Size: 1339



[View online »](#)


Cytogenetics: Xp22.31

RefSeq ORF: 558

Synonyms: TEX39B

Summary: This gene is a member of a gene family which arose through duplication on the X chromosome. The encoded protein may be localized to the nucleus as the protein contains several nuclear localization signals, and has similarity to a synaptonemal complex protein. [provided by RefSeq, Aug 2011]

Product images:



Coomassie blue staining of purified FAM9B protein (Cat# TP313736). The protein was produced from HEK293T cells transfected with FAM9B cDNA clone (Cat# [RC213736]) using MegaTran 2.0 (Cat# [TT210002]).