

Product datasheet for **TP313742L**

CTRP3 (C1QTNF3) (NM_181435) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human C1q and tumor necrosis factor related protein 3 (C1QTNF3), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213742 representing NM_181435 Red =Cloning site Green =Tags(s)

MLWRQLIYWQLLALFFLPFCLCQDEYMEVSGRTNKVVARIVQSHQQTGRSGSRREKVRERSHPKTGTVDN
NTSTDLKSLRPDELPHPEVDDLAQITTFWGGQSPQTGGLPPDCSKCCHGDYSFRGYQGPPGPPGIPGN
HGNNGNGATGHEGAKGEKGDKDLGPRGERGQHGPKEGKYPGIPPELQIAFMASLATHFSNQNSGIIF
SSVETNIGNFFDVMTGRFGAPVSGVYFFTFSMMKHEDVEEVVYLMHNGNTVFSMYSYEMKGSSTSSNH
AVLKLAKGDEWLRMGNGALHGDHQRFFSTFAGFLFETK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	35 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_852100
Locus ID:	114899



[View online »](#)

UniProt ID:	Q9BXJ4
RefSeq Size:	3789
Cytogenetics:	5p13.2
RefSeq ORF:	957
Synonyms:	C1ATNF3; CORCS; CORS; CORS-26; CORS26; CTRP3
Protein Families:	Secreted Protein

Product images:



Coomassie blue staining of purified C1QTNF3 protein (Cat# [TP313742]). The protein was produced from HEK293T cells transfected with C1QTNF3 cDNA clone (Cat# [RC213742]) using MegaTran 2.0 (Cat# [TT210002]).