

Product datasheet for TP313859

SNX29 (NM_001080530) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human sorting nexin 29 (SNX29), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213859 representing NM_001080530 Red =Cloning site Green =Tags(s)
	<p>MHSWAPLKVHLHNSDILFPVSGVGSYSPADAPLGSLNGTGPEDHVLDPGLRYSVEASSPGHGSPLSSL LPSASVPESMTISELRQATVAMMNRKDELEENRSLRNLLDGEMEHSAALRQEVDTLKRKVAEQEERQGM KVQALARENEVLKVQLKKYVGAVQMLKREGQTAEVPLNLSVDGEVTVAEQKPGEIAEELASSYERKLIEV AEMHGELIEFNRLHRLVAKEALVSQMRQELIDLRGPVPGDLSQTSQSDQSLDFEISNRALINWIPSV FLRGKAANAFHVYQVYIRIKDDEWNIYRRYTEFRSLHHKLQNKYPQVRAYNFPPKKAIGNKDAKFVEERR KQLQNYLRSVMNKVIQMVPEFAASPCKETLIQLMPFFVDITPPGEPVNSRPKAASRFKLSRGQPRETRN VEPQSGDL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	48.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001073999



[View online »](#)

Locus ID:	92017
UniProt ID:	Q8TEQ0
RefSeq Size:	1540
Cytogenetics:	16p13.13-p13.12
RefSeq ORF:	1371
Synonyms:	A-388D4.1; FLJ00143

Product images:

Coomassie blue staining of purified SNX29 protein (Cat# TP313859). The protein was produced from HEK293T cells transfected with SNX29 cDNA clone (Cat# [RC213859]) using MegaTran 2.0 (Cat# [TT210002]).