

Product datasheet for TP313946

CABP5 (NM_019855) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human calcium binding protein 5 (CABP5), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213946 protein sequence Red =Cloning site Green =Tags(s) MQFPMGPACIFLRKGIAEKQRERPLGQDEIEELREAFLEFDKDRDGFISCKDLGNLMRTMGYPTEMELI ELGQQIRMNLGGRVDFDDFVELMTPKLLAETAGMIGVQEMRDAFKEFDTNKDGEITLVELQQAMQRLLGE RLTPREISEVVREADVNGDGTVDFFEFVKMMMSR TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	19.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_062829</u>
Locus ID:	56344
UniProt ID:	<u>Q9NP86</u>
RefSeq Size:	1845


[View online »](#)


Cytogenetics: 19q13.33

RefSeq ORF: 519

Synonyms: CABP3

Summary: The product of this gene belongs to a subfamily of calcium binding proteins, which share similarity to calmodulin. Calcium binding proteins are an important component of calcium mediated cellular signal transduction. Expression of this gene is retina-specific. The mouse homolog of this protein has been shown to express in the inner nuclear layer of the retina, suggested its role in neuronal functioning. The specific function of this gene is unknown. [provided by RefSeq, Oct 2009]

Product images:



Coomassie blue staining of purified CABP5 protein (Cat# TP313946). The protein was produced from HEK293T cells transfected with CABP5 cDNA clone (Cat# [RC213946]) using MegaTran 2.0 (Cat# [TT210002]).