

Product datasheet for TP313954L

C20orf27 (NM_001039140) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 20 open reading frame 27 (C20orf27), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC213954 representing NM_001039140 Red =Cloning site Green =Tags(s) MAAANKGKCLPGVVGLAQLPVGPGRRRAIAAGNKPRVRSIRFAAGHDAEGSHSHVHFDEKLHDSVMVT Q ESDSSFLVKVGFLLKILHRYEITFTLPPVHRLSKDVREAPVPSLHLKLLSVVPVPEGYSVKCEYSAHKEGV LKEEILLACEGGTGTCVRVTVQARVMDRHHGTPMLLDGVKCVGAEEYDSEHSDWHGFD TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	21.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001034229
Locus ID:	54976
UniProt ID:	Q9GZN8



[View online »](#)

RefSeq Size: 1376

Cytogenetics: 20p13

RefSeq ORF: 597

Product images:



Coomassie blue staining of purified C20orf27 protein (Cat# [TP313954]). The protein was produced from HEK293T cells transfected with C20orf27 cDNA clone (Cat# [RC213954]) using MegaTran 2.0 (Cat# [TT210002]).