

Product datasheet for **TP314069L**

CBL (NM_005188) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Cas-Br-M (murine) ecotropic retroviral transforming sequence (CBL), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214069 representing NM_005188 Red =Cloning site Green =Tags(s)
	<p>MAGNVKKSSGAGGGSGSGSGGSLIGLMKDAFQPHHHHHHLSPPHGGTVDKKMVEKCKWLMDKVVRLC QNPKLALKNSPPYILDLLPDTYQHLRTILSRYEGKMETLGENEYFRVFMENLMKKTQKQISLFKEGKERM YEENSQPRRNLTKLSLIFSHMLAELKGIFPSGLFQGDTRITKADAAEFWRKAFGEKTIVPWKSFRQALH EVHPISSGLEAMALKSTIDLTCNDYISVFEFDIFTRLFQPWSSLLRNWNSLAVTHPGYMAFLTYDEVKAR LQKFIHKPGSYIFRLSCTRLGQWAIGYVTADGNILQTIHPNKPLFQALIDGFREGFYLPDGRNQNPDLT GLCEPTQDHIKVTQEYELYCEMGSTFQLCKICAENDKDVKIEPCGHLMCTSLTSWQESEGGQCPFCR CEIKGTEPIVDPFDPRGSGSLLRQGAEGAPSPNYDDDDDERADDTLFMMKELAGAKVERPPSPFMAPQ ASLPPVPPRLDLLPQRVCVPSSASALGTASKAASGSLHKDKPLPVPPTLRDLPPPPPPDRPYSVGAESRP QRRPLPCTPGDCPSRDKLPPVPSSRLGDSWLPPIPVPVSAPSSDPWTGRELTNRHSLPFLPSQMEP RPDVPRLGSTFSLDTSMSMNSSPLVGPECDHPKIKPSSANAIYSLAARPLPVPKLPPEQCEGEEDTEY MTPSSRPLRPLDTSQSSRACDCDQQIDSCYEAMYNISQAPSITESSTFGEGNLAAAHANTGPEESENE DDGYDVPKPPVPAVLARRTLSDISNASSFGWLSLDGDPPTTNVTEGSQVPERPPKPPRRINSERKAGSC QQGSGPAASAATASPQLSSEIENLMSQGYSDIQKALVIAQNNIEMAKNILREFVSISSPAHVAT</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	99.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005179
Locus ID:	867
UniProt ID:	P22681
RefSeq Size:	11242
Cytogenetics:	11q23.3
RefSeq ORF:	2718
Synonyms:	C-CBL; CBL2; FRA11B; NSLL; RNF55
Summary:	<p>This gene is a proto-oncogene that encodes a RING finger E3 ubiquitin ligase. The encoded protein is one of the enzymes required for targeting substrates for degradation by the proteasome. This protein mediates the transfer of ubiquitin from ubiquitin conjugating enzymes (E2) to specific substrates. This protein also contains an N-terminal phosphotyrosine binding domain that allows it to interact with numerous tyrosine-phosphorylated substrates and target them for proteasome degradation. As such it functions as a negative regulator of many signal transduction pathways. This gene has been found to be mutated or translocated in many cancers including acute myeloid leukaemia, and expansion of CGG repeats in the 5' UTR has been associated with Jacobsen syndrome. Mutations in this gene are also the cause of Noonan syndrome-like disorder. [provided by RefSeq, Jul 2016]</p>
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Chronic myeloid leukemia, Endocytosis, ErbB signaling pathway, Insulin signaling pathway, Jak-STAT signaling pathway, Pathways in cancer, T cell receptor signaling pathway, Ubiquitin mediated proteolysis

Product images:

Coomassie blue staining of purified CBL protein (Cat# [TP314069]). The protein was produced from HEK293T cells transfected with CBL cDNA clone (Cat# [RC214069]) using MegaTran 2.0 (Cat# [TT210002]).