

Product datasheet for **TP314181M**

Ankyrin brain (ANK2) (NM_020977) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ankyrin 2, neuronal (ANK2), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T



[View online »](#)

Expression cDNA Clone >RC214181 representing NM_020977
or AA Sequence: Red=Cloning site Green=Tags(s)


MMNEDAAQKSDSGEKFNQSSQRRKRPKSDSNASFLRAARAGNLDKWEYLKGGIDINTCNQNGLNALHL
AAKEGHVGLVQELLGRGSSVDSATKKGNTALHIASLAGQAEVVKVLVKEGANINAQSQNGFTPLYMAAQE
NHIDVWKYLLENGANQSTATEDGFTPLAVALQQGHNQAVAILLENDTKGKVRLPALHIAARKDDTKSAAL
LLQNDHNADVQSKMMVNRTTESGFTPLHIAAHYGNVNVATLLNLRGAADVFTARNGITPLHVASKRGNTN
MVKLLLDGRGGQIDAKTRDGLTPLHCAARSGHDQVVELLLERGAPLLARTKNGLSPLHMAAQGDHVECVKH
LLQHKAPVDDVTLDYLTALHVAHCHGHYRVTKLLLDKRAPNARALNGFTPLHIACKNRIKVMELLVKY
GASIQAITESGLTPIHVAAFMGHLNIVLLLLQNGASPDVTNIRGETALHMAARAGQVEVRCLLRNGALV
DARAREEQTPLHIASRLGKTEIVQLLLQHMAHPDAATTNGYTPLHISAREGQVDVASVLEAGAAHSLAT
KKGFTPLHVAAKYGSGLDVAKLLQRRAAADSAGKNGLTPHVAHYDNQKVALLLLEKGGASPHATAKNGY
TPLHIAAKKNQMQUIASTLLNYGAETNIVTKQGVTPHLASQEGHTDMVTLKLDKGANIHMSTKSLTSLH
LAAQEDKVNVDILTCKHGADQDAHTKLGYPPLIVACHYGNVKNMNFLLKQGANVNAKTKNGYTPLHQAAQ
QGHTHIINVLLQHGAKPNATTANGNTALAIKRLGYISVVDTLKVVTEEVTTTTTITEKHKLNPETMT
EVLVDSDEEGDDTMTGDGGEYLRPEDLKELGDDSLPSSQFLDGMNYLRYLSEGGRSDSLRSFSSDRSHTL
SHASYLRDSAVMDDSVIPSHQVSTLAKEAERNYSYRLSWGTEENLDNVALSSSPIHSGFLVFSFMDARGGA
MRGCRHNGLRIIPPRKCTAPTRVTCRLVKRHRLATMPPMVEGEGLASRLIEVGPSTGAQFLGKHLPTAP
PPLNEGESLVSRLQLGPPGKFLGPVIVEIPHFAALRGKERELVLRSENGDSWKEHFCDYTEDELNEI
LNGMDEVLDSPEDLEKKRICRIITRDFPQYFAVVSRIKQDSNLIGPEGGVLSSVTPVQVAVFPEGALTK
RIRVGLQAQPMHSELVKKILGNKATFSPIVTLEPRRRKFHKPITMTIPVPAKSSDVMLNGFGGDAPTLLR
LCSITGGTTPAQWEDITGTTPLTFVNECVSFTTNVSARFWLIDCRQIQESVTFASQVYREIICVPYMAKF
VFAKSHDPIEARLRCFCMTDDKVDKLEQQENFAEVARSRDVEVLEGGKPIYVDCFGNLVPLTKSGQHII
FSFFAFKENRPLFVKVRDTPQPCGRLSFMKEPKSTRGLVHQAIACNLNITLPIYTKESSEDQEQQEIID
MTSEKNPQDEQERIEERLAYIADHLGFSWTELARELDFTEEQIHQIRIENPNLQDQSHALLKYWLERDG
KHATDTNLVECLTKINRMDIVHLMETNTEPLQERISHSYAEIEQTITLDHSEGFVSLQEELCTAQHKQKE
EQAVSKESETCDHPPIVSEEDISVGYSTFQDGVPKTEGDSSATALFPQTHKEQVQDFSGKMQDLPEESS
LEYQQEYFVTPGTETSETQKAMIVPSSPKTPEEVSTPAEEEEKLYLQPTSSERGGSPHIIQEPEEPSEH
REESSPRKTSLVIVESADNQPETCERLDEDAAFEKGGDMPEIPPETVTEEEYIDEHGHTVWKKVTRKIIR
RYVSSEGTEKEEIMVQGMPEPVNIEEGDGYSKVIKRVWLKSDTEQSEDNNE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

- Tag:** C-Myc/DDK
- Predicted MW:** 205.6 kDa
- Concentration:** >0.05 µg/µL as determined by microplate BCA method
- Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining
- Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
- Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
- Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_066187
Locus ID:	287
UniProt ID:	Q01484
RefSeq Size:	7911
Cytogenetics:	4q25-q26
RefSeq ORF:	5616
Synonyms:	ANK-2; brank-2; CFAP87; FAP87; LQT4
Summary:	This gene encodes a member of the ankyrin family of proteins that link the integral membrane proteins to the underlying spectrin-actin cytoskeleton. Ankyrins play key roles in activities such as cell motility, activation, proliferation, contact and the maintenance of specialized membrane domains. Most ankyrins are typically composed of three structural domains: an amino-terminal domain containing multiple ankyrin repeats; a central region with a highly conserved spectrin binding domain; and a carboxy-terminal regulatory domain which is the least conserved and subject to variation. The protein encoded by this gene is required for targeting and stability of Na/Ca exchanger 1 in cardiomyocytes. Mutations in this gene cause long QT syndrome 4 and cardiac arrhythmia syndrome. Multiple transcript variants encoding different isoforms have been described. [provided by RefSeq, Dec 2011]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified ANK2 protein (Cat# [TP314181]). The protein was produced from HEK293T cells transfected with ANK2 cDNA clone (Cat# [RC214181]) using MegaTran 2.0 (Cat# [TT210002]).