

Product datasheet for TP314217L

HOXB1 (NM_002144) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human homeobox B1 (HOXB1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214217 representing NM_002144 Red =Cloning site Green =Tags(s)
	 MDYNRMNSFLEYPLCNRGPSAYSAHSAPTSFPPSSAQAVDSYASEGRYGGGLSSPAFQQNSGYPAQQPPS TLGVFPSSAPSGYAPAACSPSYGPSQYYPLGQSEGDGGYFHPSSYGAQLGGLSDGYGAGGAGPGPYPPQ HPPYGNEQTASFAPAYADLLSEDKETPCPSEPNTPTARTFDWMKVKRNPPKTAKVSEPLGSPSGLRTNF TTRQLTELEKEFHFNKYLSRARRVEIAATLELNETQVKIWFQNRMMKQKKREREGRVPPAPPGCPKEAA GDASDQSTCTSPEASPSVTS TRTRPLEQLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	32 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002135
Locus ID:	3211




[View online »](#)

UniProt ID: [P14653](#)
RefSeq Size: 1014
Cytogenetics: 17q21.32
RefSeq ORF: 903
Synonyms: HCFP3; Hox-2.9; HOX2; HOX2I

Summary: This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, located on different chromosomes, consisting of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXB genes located in a cluster on chromosome 17. [provided by RefSeq, Jul 2008]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified HOXB1 protein (Cat# [TP314217]). The protein was produced from HEK293T cells transfected with HOXB1 cDNA clone (Cat# [RC214217]) using MegaTran 2.0 (Cat# [TT210002]).