

Product datasheet for **TP314283L**

ACSM2B (NM_182617) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human acyl-CoA synthetase medium-chain family member 2B (ACSM2B), nuclear gene encoding mitochondrial protein, transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214283 representing NM_182617 Red =Cloning site Green =Tags(s)
	MHWLRKVQGLCTLWGTQMSSRTLYINSRQLVSLQWGHQEVPAKFNFASDVLDHWADMEKAGKRLPSP ALW WVNGKGKELMWNFRELENSQQAANVLSGACGLQRGDRVAVMLPRVPEWWLVILGCIRAGLIFMPGTI QM KSTDILYRLQMSKAKAIVAGDEVIQEVDTVASECPSLRIKLLVSEKSCDGLNFKLLNEASTTHHCVET GSQEASAIYFTSGTSGLPKMAEHSYSSLGLKAKMDAGWTGLQASDIMWTISDTGWILNIGSLLESWTLG ACTFVHLLPKFDPLVILKTLSSYPIKSMMGAPIVYRMLLQQDLSSYKPHLQNCLAGGESLLPETLENWR AQTGLDIREFYGQTETGLTCMVSKTMKIKPGYMGTAASCYDVQVIDDKGNVLP PGTEGDIGIRVKPIRPI GIFSGYVENPDKTAAANIRGDFWLLGDRGIKDEDGYFQFMGRADDIINSSGYRIGPSEVENALMKHPAVVE TAVISSPDPVRGEVVKAFVILASQFLSHDPEQLTKELQQHVKSVTAPYKYPRKIEFVLNLPKTVTGKIQR TKLRDKEWKMSGKARAQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	64.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_872423
Locus ID:	348158
UniProt ID:	Q68CK6
RefSeq Size:	1995
Cytogenetics:	16p12.3
RefSeq ORF:	1731
Synonyms:	ACSM2; HXMA; HYST1046
Summary:	Has medium-chain fatty acid:CoA ligase activity with broad substrate specificity (in vitro). Acts on acids from C(4) to C(11) and on the corresponding 3-hydroxy- and 2,3- or 3,4-unsaturated acids (in vitro).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified ACSM2B protein (Cat# [TP314283]). The protein was produced from HEK293T cells transfected with ACSM2B cDNA clone (Cat# [RC214283]) using MegaTran 2.0 (Cat# [TT210002]).