

## Product datasheet for TP314299

### Cardiac Troponin T (TNNT2) (NM\_001001432) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens troponin T type 2 (cardiac) (TNNT2), transcript variant 4, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214299 representing NM_001001432 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MSDIEEVVEEYEEEEQEEQEEAAEEDAEEAEAEETEETRAEEDEEEEEAKEAEDGPMEESKPKPRSFMPNLV PPKIPDGERVDFDDIHRKRMEKDLNELQALIEAHFENRKKEEELVSLKDRIERRRAERAEQQRIRNERE KERQNRLAEERARREEEENRRKAEDEARKKKALSNNMMHFGGYIQKAQTERKSGKRQTEREKKKILAERR KVLAIIDLHNLNEDQLREKAKELWQSIYNLEAEKFDLQEKFKQKYEINVLNRINDNQKVKSTRGKAKVTGR WK
	<b>TR</b> TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	33.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001001432</a>
Locus ID:	7139



[View online »](#)

UniProt ID: [P45379](#)

RefSeq Size: 1114

Cytogenetics: 1q32.1

RefSeq ORF: 846


Synonyms: CMD1D; CMH2; CMPD2; cTnT; LVNC6; RCM3; TnTC

**Summary:** The protein encoded by this gene is the tropomyosin-binding subunit of the troponin complex, which is located on the thin filament of striated muscles and regulates muscle contraction in response to alterations in intracellular calcium ion concentration. Mutations in this gene have been associated with familial hypertrophic cardiomyopathy as well as with dilated cardiomyopathy. Transcripts for this gene undergo alternative splicing that results in many tissue-specific isoforms, however, the full-length nature of some of these variants has not yet been determined. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome

**Protein Pathways:** Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM)

### Product images:



Coomassie blue staining of purified TNNT2 protein (Cat# TP314299). The protein was produced from HEK293T cells transfected with TNNT2 cDNA clone (Cat# [RC214299]) using MegaTran 2.0 (Cat# [TT210002]).