

Product datasheet for **TP314308M**

VAPA (NM_003574) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa (VAPA), transcript variant 1, 100 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC214308 representing NM_003574
Red=Cloning site **Green**=Tags(s)

MASASGAMAKHEQILVLDPPTDLKFKGPFTDVVTTNLKLRNPSDRKVCFKVKTAPRRYCVRPNNGIIDP
GSTVTVSVMLQPFDYDPNEKSKHKFMVQTI FAPPNTSDMEAVWKEAKPDELMDSKLRCVFEMPENENDKLG
ITPPGNAPTVMSSINNTVATPASYHTKDDPRGLSVLKQEKQKNDMEPSKAVPLNASKQDGPMPKPHSV
SLNDTETRKLMEECKRLQGEMMKLSEENRHLRDEGLRLRKAHSDKPGSTSTASFRDNVTSPLPSLLVVI
AAIFIGFFLGKFIL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 32.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_003565](#)


Locus ID: 9218



[View online »](#)

UniProt ID:	Q9P0L0
RefSeq Size:	6994
Cytogenetics:	18p11.22
RefSeq ORF:	882
Synonyms:	hVAP-33; VAMP-A; VAP-33; VAP-A; VAP33
Summary:	The protein encoded by this gene is a type IV membrane protein. It is present in the plasma membrane and intracellular vesicles. It may also be associated with the cytoskeleton. This protein may function in vesicle trafficking, membrane fusion, protein complex assembly and cell motility. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Jul 2008]
Protein Families:	Transmembrane
Protein Pathways:	Tight junction

Product images:



Coomassie blue staining of purified VAPA protein (Cat# [TP314308]). The protein was produced from HEK293T cells transfected with VAPA cDNA clone (Cat# [RC214308]) using MegaTran 2.0 (Cat# [TT210002]).