

Product datasheet for TP314311L

HOMER2 (NM_199331) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens homer homolog 2 (Drosophila) (HOMER2), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214311 representing NM_199331 Red =Cloning site Green =Tags(s)
	MGEQPIFTTRAHVVFQIDPNTKKNWMPASKQAVTVSYFYDVTRNSYRIISVDGAKVIINSTITPNMTFTKT SQKFGQWADSRANTVFGGLGFSSEQQLTKFAEKFEVKEAAKIAKDKTQEKIETSSNHSQASSVNGTDDEK ASHAGPANTHLKSENDKLIKIALTQSAANVKNEINREKEKNTQLKRRIEELEAELREKETELKDLRKQSEI IPQLMSECEYVSEKLEAAERDNQNLEDKVRSLKTDIEESKYRQRHLKVELKSFLEVLDGKIDDLHDFRRG LSKLGTDN
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	32.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_955363
Locus ID:	9455



[View online »](#)

UniProt ID:	Q9NSB8
RefSeq Size:	1807
Cytogenetics:	15q25.2
RefSeq ORF:	864
Synonyms:	ACPD; CPD; HOMER-2; HOMER2A; HOMER2B; Vesl-2
Summary:	This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function. The encoded protein is a postsynaptic density scaffolding protein. Alternative splicing results in multiple transcript variants. Two related pseudogenes have been identified on chromosome 14. [provided by RefSeq, Jun 2011]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified HOMER2 protein (Cat# [TP314311]). The protein was produced from HEK293T cells transfected with HOMER2 cDNA clone (Cat# [RC214311]) using MegaTran 2.0 (Cat# [TT210002]).